STERILE M.U.S.T.[®] and M.U.S.T.[®] MC MODULAR TAP

CAUTION

For a U.S. specific package insert, please go to www.medacta.com or call toll free 800-901-7836.
 All Medacta M.U.S.T. and M.U.S.T. MC Sterile Instruments are supplied in single-use packages. The sterilization method is indicated on the label

SYMBOLS / SIMBOLI / PICTOGRAMMES / PIKTOGRAMME / PICTOGRAMAS / PICTOGRAMAS / PICTOGRAMMES

	Do not reuse / Non riutilizzare / Niet hergebruiken / Nicht wiederverwenden / No reutilizable / Não reutilizar / Ne pas réutiliser
	Do not resterilize / Non risterilizzare / Niet reesteriliseren / Bitte nicht re-sterilisieren / No reesterilize / Não reesterilizar / Ne pas restériliser
	Caution, read the accompanying documents / Attenzione, consultare i documenti di accompagnamento / Opgelet, raadpleeg de bijgeleverde documenten / Achtung, Begleitdokumente beachten / Atención, consulte los documentos que se acompañan / Atenção, consultar os documentos de acompanhamento / Attention, consulter les documents d'accompagnement
	Consult instructions for use / Consultare le istruzioni di utilizzo / Gebruiksaanwijzing / Bitte lesen Sie die Gebrauchsanweisung / Consulte las Instrucciones antes de su uso / Consulte as instruções de uso / Consultez les instructions d'utilisation
	Do not expose to sunlight / Non esporre alla luce del sole / Niet blootstellen aan zonlicht / Vor Sonnenlicht schützen / No exponer a la luz solar / Não expor à luz solar / Ne pas exposer à la lumière du soleil
	Store in a dry place / Conservare in luogo asciutto / Droog bewaren / Trocken aufbewahren / Conservar en lugar seco / Conservar em local seco / Conserver au sec
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Niet gebruiken als verpakking beschadigd is / Bitte nicht verwenden, falls Verpackung beschädigt ist / No utilize si el envoltorio está dañado / Não utilizar se a embalagem estiver violada / Ne pas utiliser si l'emballage est abîmé
	Use by / Utilizzare entro il / Te gebruiken voor / Zu verwenden bis / Utilizar antes de / A utilizar antes de / A utiliser avant
	Lot number / Numero di lotto / Nummer van de partij / Waren-Lot-Nummer / Número de lote / Número de lote / Numéro de lot
	Reference number / Riferimento commerciale / Handelsreferentie / Referenznummer / Numero de Referencia / Referência comercial / Référence commerciale
	Sterilized with ethylene oxide / Sterilizzato con ossido di etilene / Gesteriliseerd met ethyleenoxide / Ethylenoxid-sterilisiert / Esterilizado con óxido de etileno / Esterilizado com óxido de etileno / Stérilisé à l'oxyde d'éthylène
	Sterilized by irradiation / Sterilizzato mediante irradiazione / Gesteriliseerd door doorstraling / Sterilisiert durch Bestrahlung / Esterilizado por irradiación / Esterilizado por irradiação / Stérilisé par irradiation
	Manufacturer / Produttore / Fabrikant / Hersteller / Fabricante / Fabricante / Fabricant
	Double sterile barrier system / Sistema a Doppia Barriera Sterile / Dubbele steriele verpakking / Doppeltes steriles Barriersystem / Sistema de doble barrera estéril / Sistema de Dupla Barreira Estéril / Système à barrière stérile double / Σύστημα διπλού στείρου εμποδίου

Index

ENGLISH: M.U.S.T. AND M.U.S.T MC MODULAR TAP - INSTRUCTIONS FOR USE	3
ITALIANO: M.U.S.T. E M.U.S.T. MC MODULAR TAP - ISTRUZIONI PER L'USO	5
NEDERLANDS: M.U.S.T. EN M.U.S.T. MC MODULAR TAP - GEBRUIKSINSTRUCTIES	7
DEUTSCH: M.U.S.T. UND M.U.S.T. MC MODULAR TAP - GEBRAUCHSANWEISUNGEN	9
ESPAÑOL: M.U.S.T. Y M.U.S.T. MC MODULAR TAP - INSTRUCCIONES DE USO	11
PORTUGUÊS: M.U.S.T. E M.U.S.T. MC MODULAR TAP - INSTRUÇÕES DE UTILIZAÇÃO	13
FRANÇAIS : M.U.S.T. ET M.U.S.T. MC MODLAR TAP - INSTRUCTIONS D'UTILISATION.....	15

ENGLISH: M.U.S.T. AND M.U.S.T MC MODULAR TAP - INSTRUCTIONS FOR USE

Important notice: the device(s) can be prescribed and implanted only by a doctor legally authorized to perform this type of surgery.

1. GENERAL

Use of the instrumentation requires knowledge of the anatomy and pathology, biomechanics and surgical corrections of the spine. The instrumentation can be used only by a qualified surgeon who practices with an awareness of current advances in science and surgical techniques. The user should ensure that the instruments are intact and in good working order before use. The user should also take all necessary precautions to avoid any accidents (gloves, protective glasses...). No undesirable side effects are known if these Instructions for use are respected. Before any surgery, the surgeon must be familiar with the sales product literature and operative technique and must carefully read these Instructions for use. The instrumentation should be used only for its intended purpose as indicated in the surgical techniques. The use of some motor-driven instruments (drill bits, taps, etc.) may cause a temperature rise between these instruments and the bone. It is advisable to spray these instruments with physiological saline solution during use. Under no circumstances should an item of instrumentation be implanted.

2. PRODUCT DESCRIPTION

The instrumentation is made from materials appropriate for the manufacture of surgical instruments. These materials are not intended to stay in permanent contact with the patient. The instrumentation is supplied sterile and it is intended to be for single use. The instrumentation manufactured by Medacta International meets the mechanical and functional characteristics of the operative technique and these instructions. Before surgery, the user should refer to the surgical technique and other labeling or contact the representative of the company for more details on how to use the instrumentation.

3. INTENDED USE

The modular taps are part of the M.U.S.T. and M.U.S.T MC pedicle screw system, which is intended to provide stabilization and to promote bone fusion during the normal healing process following the surgical correction of disorders of the spine. Specifically, the modular taps are intended to prepare the pedicle screw canal intra-operatively.

4. INDICATIONS

Modular taps are intended for use with the M.U.S.T. and M.U.S.T MC pedicle screw system and according to its approved indications for use.

5. CONTRAINDICATIONS

Contraindications for using modular taps are the same as contraindications for a spinal fusion with pedicle screws. Please refer to the M.U.S.T. and M.U.S.T MC surgical technique for a comprehensive list of contraindications.

6. WARNINGS

Check the expiration date prior to use. Care should be taken to ensure that the correct size (diameter and length) is being used. At any time during the surgical procedure, the surgeon can verify and confirm that the modular taps are positioned correctly on the pedicles by fluoroscopy. Visually inspect the instruments after use in order to verify they did not experience any mechanical damage which may cause release of particles in the human body. Any malfunctioning instruments should be immediately returned to Medacta. The type of malfunction should also be reported.

7. HANDLING

The instrumentation should be handled with care. The surgical instruments can be damaged by inappropriate handling; visually inspect the instrument and check for damage prior to use.

8. PACKAGING

The modular taps are supplied sterile, in single-use individual packages. The sterilization method is indicated on the label. The expiration date and package integrity must be checked in order to ensure that the sterility of the contents has not been compromised. If the package is damaged, do not use the component. Do not resterilize.

9. STORAGE

The packages must be stored in a cool, dry place, not exposed to sunlight.

10. FURTHER INFORMATION

Recommended directions for use of this system (surgical operative techniques) are available at no charge upon request. If further information is needed or required, please contact Medacta.

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