

# U.S.T. MINI OPEN

MINIMALLY INVASIVE RETRACTOR SYSTEM

M.I.S. MAKE IT SMART IN ONE SYSTEM



## Surgical Technique

Hip

Knee

Spine

Navigation




## C O N T E N T S

<b>1</b>	<b>INTRODUCTION</b>	<b>4</b>
<b>2</b>	<b>SURGICAL TECHNIQUE</b>	<b>5</b>
<b>3</b>	<b>SURGICAL STEPS FOR MINI OPEN RETRACTION</b>	<b>5</b>
3.1	Pedicle Preparation	5
3.2	Dilator Insertion	5
3.3	Retractor Placement	6
3.4	Retractor System Lock	7
3.5	Frame Expansion	8
3.6	Enhanced Retraction	8
3.7	Light System	9
<b>4</b>	<b>SURGICAL STEPS FOR PEDICLE SCREW INSERTION</b>	<b>10</b>
4.1	Screw Insertion	10
4.2	Parallel Distraction	10
<b>5</b>	<b>SURGICAL STEPS FOR IMPLANT-BLADE HYBRID TECHNIQUE</b>	<b>12</b>
5.1	Pedicle Preparation	12
5.2	Pedicle Screw and Blade Placement	12
5.3	Percutaneous Tower Removal	13
<b>6</b>	<b>SURGICAL STEPS FOR INTERBODY FUSION DEVICE INSERTION</b>	<b>14</b>
<b>7</b>	<b>SURGICAL STEPS FOR ROD PLACEMENT AND ADJUSTMENTS</b>	<b>15</b>
7.1	Rod Insertion	15
7.2	Compression and Distraction	15

## 1 INTRODUCTION

The Medacta MINI OPEN Retractor allows for decompression, fusion and fixation through a minimally invasive approach, however providing the surgeon the same versatility of an open procedure.


The Medacta MINI OPEN Retractor represents a gentle interface to the soft tissues, potentially allowing to reduce blood loss and scarring, improve the patient's recovery time and shorten the hospital stay.



### Features

#### Cranial / Caudal blades

- Radiolucent (Aluminium)
- Left / Right Configuration
- Anatomical blades. With a specific design matching the spine anatomy. 60 – 100 mm Length (10 mm increment)
- MIS Blades. Straight design. 50-120mm (10 mm increment)
- Illuminated anatomical blades. To be used in combination with a light source and cables with a specific design matching the spine anatomy. 30 – 100 mm Length (10 mm increment, right and left configuration)
- Illuminated MIS blades, straight design. To be used in combination with a light source and cables. 30 – 120 mm Length (10 mm increment)



#### Medial / Lateral blades

- Radiolucent (Aluminium)
- Convex and Straight Design
- 50 – 120mm Length (10 mm increment)
- 15mm and 20mm wide configuration



#### Cranial / Caudal Frame

- Couple with the medial / lateral blades for optimal distraction
- Allows for blades 30° tilting excursion for additional access

#### Medial / Lateral Frame

- Fits under the cranial / caudal blades for optimal distraction

#### MIS Medial Arms

- Couple with the Cranial / Caudal frame for aided distraction / retraction
- Allows for blades 30° tilting excursion for additional access
- Two designs available: Ratcheting and Pivoting



#### Distraction Handle

- Medial and Lateral
- Allows medial / lateral independent distraction

### Light cables

- Bifurcated surgical illuminator. Designed to be mounted on the illuminated blades, packaged sterile
- Cable extension, to connect the light source to the bifurcated surgical illuminator
- Adaptors, to connect the light source to the cable extension




## 2 SURGICAL TECHNIQUE

Perform a posterior incision through the skin and fascia approximately 2 – 4 cm from the midline. Allow for muscle splitting within the multifidus and longissimus cleavage plane.

## 3 SURGICAL STEPS FOR MINI OPEN RETRACTION


### 3.1 Pedicle Preparation

Target the pedicle and perforate the outer cortex with the Cannulated Awl.



#### ! WARNING

Confirm the pedicle anatomical positioning with radiographic imaging.



Remove the inner pin from the Cannulated Awl and insert the K-wire, carefully advance it through the incision.


#### ! WARNING

Control the k-wire position under radiographic imaging and make sure it does not slip off during the procedure.

### 3.2 Dilator Insertion

#### Multi-Step Dilation


Insert the Dilator size 1 – 8mm over the k-wire. Continue dilation sliding over the Dilator size 2 – 16mm and finally place the Dilator size 3 – 22mm over each other.



**NOTICE:** The markings on the Dilators indicate the appropriate length of the blades. The Retractor has a diameter of 22 mm matching with the Dilator size 3.

#### Single Step Dilation - Anatomical

Alternatively, place the Anatomical Dilator aligned in the direction of muscle fibers; Rotate the Dilator thus allowing the muscle splitting.



**NOTICE:** The markings on the dilators indicates the appropriate length of the blades.

### 3.3 Retractor Placement

Select the retraction blades in accordance with the markings on the final dilator (Fig. 1). Slide the blades over the Dilator and place it through the incision until in a stable position (Fig.2).

Alternatively, attach the Cranial/Caudal blades to the Retractor Frame and insert it over the Dilator and into the incision.




Fig.1

To assemble the Frame onto the Cranial / Caudal Blades, the Frame must be initially set to the "unlocked" configuration with dedicated T-Handle (Fig. 2).




Fig.2

**NOTICE:** When set to the "unlocked" position the two arms are free to slide apart to facilitate engagement with the blades.

After appropriate engagement with the Blades turn the Frame to the "lock" position with the dedicated T-Handle.




Fig.3

**NOTICE:** When set to the "lock" position the two arms are engaged to the ratchet mechanism and the knob can be used for a fine-tuned longitudinal sliding.





Fig.4

**OPTION**  
 If the enhanced frame is used, the arm angulation can be locked with the THandle (Fig.4)



### 3.4 Retractor System Lock

When the Retractor System is assembled and placed over the bony anatomy, lock it in place with the Instrument Holder.


To assemble the Instrument Holder to the Retractor, align the Arm Connector to the recess of the Frame and tighten the knob with the THandle. Fix the central knob when desired angulation of the Flexible Arm is achieved.



**NOTICE:** Due to the specific coupling of the Arm Connector and the Frame, the support system grants optimal stability of the Retractor and can be assembled in three possible positions (A, B, C) to provide great flexibility at the Retracting system.



**OPTION**  
 The Arm Connector M4 can be used to fix the frame on the opposite side, coupled with the rack.




Once the frame is inserted and stable in position, the dilators and the k-wire are removed establishing a corridor for the surgery.

**OPTION**


The k-wire can be maintained for subsequent guided insertion of a cannulated screw

**3.5 Frame Expansion**

To expand the Retractor turn the paddle in a counter-clockwise direction.



Additional retraction can be reached by tilting the blades with the T-handle. Blades can be angulated up to 30° providing an aid in the visualisation of the intraoperative field (see figure below).




**3.6 Enhanced Retraction**

Medial / Lateral Blades can be used for further retraction.

Choose the length of the Medial / Lateral Blades consistently to the length of the already placed Cranial/Caudal Blades.

**Medial Retractor Blade**

Couple the preferred MIS Medial Arm with the Medial and match the assembly onto the Cranial/Caudal Frame. Retract and tilt the Blade with the T-handle, if necessary



When the Ratcheting Medial Arm is used (Fig.5), blade retraction can be performed through the proximal pin (A), while blade tilt can be adjusted by acting on the distal pin (B).






Fig.5




The Pivoting Medial Arm can also be used (Fig.6); to fix it over the frame the middle lateral pin can be regulated (C). Arm retraction and agulation can be freely adjusted and then fixed with the upper proximal pin (D), while blade tilt can be regulated by means of the distal pin (E).



### Medial And Lateral Retractor Blades

Place the Medial-Lateral Frame into the Cranial / Caudal Frame and into the incision; expand the Medial/Lateral Frame to the desired retraction length.



Alternatively, the Blade Holder can be used in combination with the Medial Arm to perform medial/lateral retraction.


The Blade Holder can also be used in combination with the Cranial/Caudal Frame to perform lateral retraction only.

### 3.7 Light System

The Mini Open set includes the illuminated blades that give the surgeon the option of an illuminated view of the situation.

The light source (not provided by the system) is attached to an extension cable (reusable) and then bifurcated cables (sterile, single use) are linked and mounted on the specifically designed illuminated blades.

The adaptors, compatible with the majority of light sources on the market, connect the light source to the Medacta light system.



## 4 SURGICAL STEPS FOR PEDICLE SCREW INSERTION

### 4.1 Screw Insertion

Using the surgical access generated with the Retractor proceed to insert the pedicle screws.  
To implant the screws in the pedicles, please follow the procedure described in the Medacta M.U.S.T. Pedicle Screw System surgical technique (ref. 99.46.12).




Fig.8

#### OPTION

The Blade-Implant Connector insertion can be guided with the k-wire when cannulated screws are used.




Fig.9

### 4.2 Parallel Distraction

The M.U.S.T. Pedicle Screw System provides the possibility to perform parallel distraction using the Retractor system; this allows the proper distraction to facilitate the interbody fusion device implantation.




Fig.7

Using the Blade-Implant Connector Holder or the Temporary Set Screwdriver, slide the Blade-Implant Connector within the blade rail to reach the tulip (Fig.7 & Fig.8).


#### WARNING

During insertion of the Blade-Implant Connector, verify the proper engagement of the dove-tail end into the rail (Fig.9).

**NOTICE:** During this phase it can be useful to switch the Frame in the "unlocked" configuration. The "lock" configuration must be recovered as soon as this step is accomplished.

Lock the connector to the Screw head using the Temporary Set Screw.

The Temporary Set screw can be inserted with the Cannulated Setscrewdriver, guided by the k-wire, or with a standard or enhanced Setscrewdriver.



Perform desired distraction by means of the frame ratchet.



Lordotic distraction can be reached by tilting the blades into a desired angle (see figure below).



## 5 SURGICAL STEPS FOR IMPLANT-BLADE HYBRID TECHNIQUE

As an alternative, it is also possible to insert the Pedicle Screw first, by a percutaneous approach, having the screw already coupled with a dedicated connector for subsequent blade and frame coupling. For this technique, cannulated screws and percutaneous tubes are used. Please refer also to 99.perc46.12. The "hybrid" technique has been developed in order to potentially save surgical time by reducing the surgical steps in comparison with the standard technique.


### 5.1 Pedicle Preparation

Prepare the pedicle cannulated Awl as described in paragraph 3.1. Remove the Awl and leave the K-wire in place. Proceed with the assembly of the enhanced cannulated pedicle screw with the Percutaneous Tower (see M.U.S.T. Percutaneous Minimal Invasive System Surgical Technique 99.perc46.12).




### 5.2 Pedicle Screw and Blade Placement


Insert the pedicle screw on the K-wire using a cannulated handle from the M.U.S.T. set.




Insert the hybrid connector 03.52.10.0217 into the Tower slot and the hybrid screw driver 03.52.10.0216 into the upper axial hole of the Tower itself.



Remove the screw driver keeping the hybrid connector in place. Insert the temporary set screw 03.52.10.0212 with the aid of the cannulated temporary set screw driver and fix, without tightening, the connector.




Slide the blade through the connector. The Percutaneous Tower acts as a dilator.




**NOTICE:** The Percutaneous Tower is provided with a measurement scale to indicate the length of the blade to be used.

### 5.3 Percutaneous Tower Removal


Remove the Percutaneous Tower by pushing the release button. After that the Kwire can be pulled out.



Repeat the previous steps in order to place the second pedicle screw along with the blade.



Proceed with the frame assembly. To fasten the engagement, the frame should be in the "unlock" position (see paragraph 3.3).



## 6 SURGICAL STEPS FOR INTERBODY FUSION DEVICE INSERTION

Using the surgical access provided with the Retractor, proceed with the interbody fusion device implantation (Fig. 10).

Perform the discectomy and the endplate preparation with dedicated instruments according to the usual fashion.

If you are planning to implant the MectalIF cages, please follow the procedure described in the dedicated surgical techniques (ref. 99.44.12 & 99.44TLIF.12)

The MectalIF Oblique and Posterior as well as the MectalIF Transforaminal Interbody Devices can be implanted in combination with the use of the M.U.S.T. Mini Open system through a minimally invasive procedure.

In particular, the MectalIF Oblique Interbody Device, which is implanted at a 30° angle in the sagittal plane, results as an optimal interface with the current retractor / distractor system (Fig. 10).



Fig. 10

The MectalIF Oblique and Posterior cages can also be implanted by using the MectalIF Posterior Offset Handle and the MectalIF Oblique Offset Handle, respectively.

The MectalIF offset handles can be coupled with either the Inner Shaft or the Extended Inner Shaft, keeping optimal visibility over the surgical space (Fig. 11).

The mark "MEDIAL" on the shaft indicates proper alignment of the instrument with respect to the patient.

Attach the implant perpendicular to the Inner Shaft / Handle assembly, screw the thread of the Inner Shaft to the threaded hole and secure it firmly.

Ensure that the orientation of the implant is correct, as shown in Fig. 12 for Oblique Implants or Fig. 13 for Posterior Implants.



Fig. 11



Fig. 12




Fig. 13


## 7 SURGICAL STEPS FOR ROD PLACEMENT AND ADJUSTMENTS

### 7.1 Rod Insertion

Using the Temporary Set Screwdriver, remove the Temporary Setscrews from the tulip and then proceed to remove the Blade-Implant Connector, with the Blade-Implant Connector Holder.

Position the rod into the selected pedicle screw heads (Fig. 14). The use of the Rod Insertion Forceps may help to easily insert the rod in the space generated between the Cranial/Caudal Blades.

Perform the Rod reduction according to the usual Techniques described in the Medacta M.U.S.T. Pedicle Screw System surgical technique (ref. 99.46.12).



### 7.2 Compression and Distraction

Further compression / distraction can be performed at this time. In either maneuver, the set screw on one side of the motion segment should be tightened while the other must be in place but untightened. The untightened screw must be coupled with the Countertorque / Set Screwdriver system. Compression occurs by acting against the tightened screw with the Compressor Lever (Fig. 15).



Similarly, distraction is reached pushing the Distraction Lever against the tightened screw (Fig. 16).




Once desired compression or distraction has been achieved, final tightening has to be performed.

For the M.U.S.T. pedicle screws, please follow the same procedure previously described in the dedicated surgical technique of the Medacta M.U.S.T. screws (99.46.12).



Finally, remove the MIS frame and the retraction blades to complete the implant procedure.









Part numbers subject to change.

---

## NOTE FOR STERILISATION

The instrumentation is not sterile upon delivery. It must be cleaned before use and sterilised in an autoclave respecting the regulations of the country, EU directives where applicable and following the instructions for use of the autoclave manufacturer.

For detailed instructions please refer to the document "Recommendations for cleaning decontamination and sterilisation of Medacta International orthopedic devices" available at [www.medacta.com](http://www.medacta.com).



---

**REDEFINING BETTER  
IN ORTHOPAEDICS  
AND NEUROSURGERY**

---

MEDACTA.COM



**Medacta International SA**  
Strada Regina - 6874 Castel San Pietro - Switzerland  
Phone +41 91 696 60 60 - Fax +41 91 696 60 66  
info@medacta.ch

Find your local dealer at: [medacta.com/locations](http://medacta.com/locations)

All trademarks and registered trademarks are the property of their respective owners.

M.U.S.T. Mini Open  
Surgical Technique

ref: 99.mini46.12  
rev. 04

Last update: July 2017  
**CE 0476**