Strada Regina
CH - 6874 Castel San Pietro - Switzerland
Phone +41 91 696 60 60 - Fax +41 91 696 60 66
info@medacta.ch - www.medacta.com

Manufactured by:
Medacta International SA
CH-6874 Castel San Pietro - Switzerland

M-ARS ACL – EXTRACORTICAL FEMORAL BUTTON

CAUTION

For U.S. specific instructions for use, please go to www.medacta.com or call toll free 800-901-7836.

SYMBOLS / SIMBOLI / PICTOGRAMMES / PIKTOGRAMME / PICTOGRAMAS / PICTOGRAMAS / PICTOGRAMMES / ΕΙΚΟΝΟΣΥΜΒΟΛΑ

	Do not reuse / Non riutilizzare / Niet hergebruiken / Nicht wiederverwenden / No reutilizable / Não reutilizar / Ne pas réutiliser / Μην επαναχρησιμοποιείτε
	Do not resterilize / Non risterilizzare / Niet resteriliseren / Bitte nicht re-sterilisieren / No resterilize / Não resterilizar / Ne pas restériliser / Να μην αποστειρωθεί ξανά
	Caution, read the accompanying documents / Attenzione, consultare i documenti di accompagnamento / Opgelet, raadpleeg de bijgeleverde documenten / Achtung, Begleitdokumente beachten / Atención, consulte los documentos que se acompañan / Atenção, consultar os documentos de acompanhamento / Attention, consulter les documents d'accompagnement / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα
	Consult instructions for use / Consultare le istruzioni di utilizzo / Gebruiksaanwijzing / Bitte lesen Sie die Gebrauchsanweisung / Consulte las instrucciones antes de su uso / Consulte as instruções de uso / Consultez les instructions d'utilisation / Οδηγίες χρήσης
	Do not expose to sunlight / Non esporre alla luce del sole / Niet blootstellen aan zonlicht / Vor Sonnenlicht schützen / No exponer a la luz solar / Não expor à luz solar / Ne pas exposer à la lumière du soleil / Μην εκθέτετε σε ηλιακό φως
	Store in a dry place / Conservare in luogo asciutto / Droog bewaren / Trocken aufbewahren / Conservar en lugar seco / Conservar em local seco / Conserver au sec / Φυλάσσετε σε ξηρό χώρο
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Niet gebruiken als verpakking beschadigd is / Bitte nicht verwenden, falls Verpackung beschädigt ist / No utilize si el envoltorio está dañado / Não utilizar se a embalagem estiver violada / Ne pas utiliser si l'emballage est abîmé / Να μην χρησιμοποιηθεί εάν η συσκευασία είναι κατεστραμμένη
	Use by / Utilizzare entro il / Te gebruiken voor / Zu verwenden bis / Utilizar antes de / A utilizar antes de / A utiliser avant / Ημερομηνία λήξης
	Lot number / Numero di lotto / Nummer van de partij / Waren-Lot-Nummer / Número de lote / Número de lote / Numéro de lot / Αριθμός παρτίδας
	Reference number / Riferimento commerciale / Handelsreferentie / Referenznummer / Numero de Referencia / Referência comercial / Référence commerciale / Νούμερο παραπομπής
	Sterilized by irradiation / Sterilizzato mediante irradiazione / Gesteriliseerd door doorstraling / Sterilisiert durch Bestrahlung / Esterilizado por irradiación / Esterilizado por irradiação / Stérilisé par irradiation / Αποστειρωμένο με ακτινοβολία
	MR Conditional / A compatibilità condizionata con la risonanza magnetica / MR-conditioneel / Bedingt MR tauglich / Compatibles con RM bajo ciertas condiciones / Qualificado para ressonância magnética / RM conditionnelle / Μαγνητικός συντονισμός υπό προϋποθέσεις

Index

ENGLISH: EXTRACORTICAL FEMORAL BUTTON - INSTRUCTIONS FOR USE.....	3
ITALIANO: EXTRACORTICAL FEMORAL BUTTON - ISTRUZIONI PER L'USO.....	6
NEDERLANDS: EXTRACORTICAL FEMORAL BUTTON - GEBRUIKSINSTRUCTIES	9
DEUTSCH: EXTRACORTICAL FEMORAL BUTTON - GEBRAUCHSANWEISUNGEN	12
ESPAÑOL: EXTRACORTICAL FEMORAL BUTTON - INSTRUCCIONES DE USO	15
PORTUGUÊS: EXTRACORTICAL FEMORAL BUTTON - INSTRUÇÕES DE UTILIZAÇÃO.....	18
FRANÇAIS: EXTRACORTICAL FEMORAL BUTTON - INSTRUCTIONS D'UTILISATION.....	21
ΕΛΛΗΝΙΚΑ: M-ARS ACL – EXTRACORTICAL FEMORAL BUTTON - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ.....	24

ENGLISH: EXTRACORTICAL FEMORAL BUTTON - INSTRUCTIONS FOR USE

Important notice: the device(s) can be prescribed and implanted only by a medical doctor legally authorized to perform this type of surgery.

1. GENERAL

Before any surgery, the surgeon must be familiar with the product literature and operative technique and must read carefully these instructions for use. Patient selection is as important as implant placement or positioning. Unsuitable functional requirements may contribute to reduce the implant life. The warnings must be heeded, and the instructions for use must be strictly followed.

2. PRODUCT DESCRIPTION

These instructions for use are intended for all products described here below.

Medacta extracortical fixation button provides the orthopedic surgeon a means of accurate extracortical fixation in reconstruction surgery. The extracortical fixation button allows for endoscopic ligament reconstruction without an ancillary lateral incision.

3. INTENDED USE / INDICATIONS

The extracortical femoral button is used for reconstructive treatment of ligament ruptures as well as for femoral fixation of an implanted cruciate ligament replacement.

Reconstructive therapy of ruptures to the anterior and posterior cruciate ligament by means of autologous grafts.

4. CONTRAINDICATIONS

- Osteoporosis and osteomalacia
- Degenerative osteopathies
- Adiposity or patient overweight, leading to excessive strain on the fixation button
- Osteomata in the area in which the fixation button is to be placed
- Deformities of the bone, or general conditions of the bones which exclude implantation of a fixation button in the opinion of a physician
- Systemic diseases and metabolic disorders that may compromise the output of the surgery

5. WARNINGS AND PRECAUTIONS

The success of the operation depends on compliance with the operative technique supplied, and the proper use of the dedicated instruments specially designed for that range of implants. Malpositioning may reduce implant longevity and lead to early implant failure.

One Medacta label for each device must be properly applied to the Medacta Implant Card provided, and delivered to the patient.

MRI Compatibility

Non-clinical testing has demonstrated titanium implants are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions. Static magnetic field of 1.5 Tesla and 3 Tesla, with :

- Maximum spatial field gradient of 12,800 G/cm (128 T/m)
- Maximum force product of 231,000,000 G²/cm (231 T²/m)

Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode). Under the scan conditions defined above, titanium implants expected to produce a maximum temperature rise of less than:

- 1.5°C (2 W/kg, 1.5 Tesla) RF-related temperature increase with a background temperature increase of ≈1.3°C (2 W/kg, 1.5 Tesla)
- 3.0°C (2 W/kg, 3 Tesla) RF-related temperature increase with a background temperature increase of ≈2.3°C (2 W/kg, 3 Tesla)

after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 20.5 mm from the titanium implant when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

6. MEDACTA INTERNATIONAL IMPLANTS

Medacta International is not responsible for the use of its implant components in combination with a component from another manufacturer (unless otherwise specified by Medacta International in the surgical technique), therefore we advise against such a use. The components of a M-ARS ACL Extracortical Femoral Button should never be reimplanted. While an implant may appear undamaged, microscopic imperfections may occur and cause implant failure. The operating surgeon must be aware that even a very small superficial damage, caused for instance by a sharp tool or electrocautery, can have an influence on the endurance of the device and can lead to fracture. Medacta M-ARS ACL products are non-pyrogenic.

7. RISK FACTORS

The following risk factors, individually or together, may result in poor clinical outcomes:

- Inadequate bone quality, (e.g. osteoporosis, previous cruciate ligament operation)
- Systemic diseases or metabolic disorders
- History of infections or recurrent falls
- Drug dependence and abuse of alcohol and medicaments
- Mental incapacity of patient to understand the instructions of the physician and to comply with them
- Local bone tumors

8. INSTRUCTIONS FOR USE

PREOPERATIVE PHASE

The surgeon should verify possible patient physical limitations and mental deficiencies and he should also discuss with the patient all the details of the procedure and implant. The discussion should consider the limitations of the procedure and the constraints imposed by the selected implant. The factors which could limit the performance and stability of the implant, e.g. level of activity, should be set out to improve the patient's chances to avoid complications. The necessity to follow the postoperative instructions given by the surgeon should be fully understood by the patient. A stock of sterile implants of suitable sizes must be available and checked by the operator before surgery.

9. HANDLING

To avoid scratching or damaging the implants, these should be handled with the utmost care by qualified personnel and in an environment where conditions of hygiene are controlled. The implants should be kept in their undamaged packages.

10. SURGICAL TECHNIQUE

The surgeon should be fully familiar with the surgical technique. Supplementary information about the surgical techniques (brochure and video) and products are available on request. Careful preoperative planning, documented by X-rays, is essential.

11. POSTOPERATIVE CARE AND FOLLOW-UP

The surgeon should caution the patients to control their level of activity and avoid excessive loads on the operated limb. Moreover the surgeon should make the patients aware of the precautions to be taken in terms of exercises, treatments and limitations on activities, any limitations reported on the label, as well as exposure to magnetic fields. The patient must be told that implants can affect the results of computer tomography (CT) or magnetic resonance imaging (MRI) scans.

Periodic follow-up and X-rays are recommended to make comparisons with the immediate postoperative condition and anticipate implant related complications. Excessive physical activity, and operated limb traumas may cause early failure of the implant. If the case occurs, it is necessary to place the patient under supervision, evaluate the possible progression of the deterioration, and weigh the benefit of early revision.

12. ADVERSE EFFECTS AND COMPLICATIONS

Adverse effects that can occur in reconstructive treatment of ligament ruptures include:

- Infection, both deep and superficial
- Allergies, mild inflammatory and foreign body reactions to implant material

Some adverse effects can ultimately lead to death.

General complications include:

- Venous thrombosis with/without pulmonary embolism;
- Cardiovascular or pulmonary disturbances;
- Haematomas;
- Systemic allergic reactions;
- Systemic pain.

13. PACKAGING

The components of all M-ARS ACL Extracortical Femoral Button are supplied in single-use packages. For components delivered sterile, the sterilization method is indicated on the label. The expiration date must be checked on the label as well as the package integrity to ensure that sterility of the contents has not been compromised. If the package is damaged or has been previously opened, do not use the component. Do not resterilize.

14. STORAGE

The packages must be stored in a cool, dry place, away from light.

15. TRADEMARKS

Medacta® is registered trademark of Medacta International SA, Castel San Pietro, Switzerland.

The reference text is the English text.

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