






<span></span>	ref:75.09.046 rev.04
<b>Medacta</b>	<b>CE 2797</b>
International <span></span>	
<b>Strada Regina</b> <b>CH - 6874 Castel San Pietro - Switzerland</b> <b>Phone<span> </span>: +4 1 91 696 60 60 Fax +41 91 696 60 66</b> <b>Info@medacta.ch – www.medacta.com</b>	Manufactured by: Medacta International SA CH-6874 Castel San Pietro - Switzerland

<b>MECTASCREW PEEK</b>
The components of all MectaScrew PEEK interference screw are supplied in single-use packages. For components delivered sterile, the sterilization method is indicated on the label. The expiration date must be checked on the label as well as the package integrity to ensure that sterility of the contents has not been compromised. If the package is damaged or has been previously opened, do not use the component. Do not resterilize.
<b>CAUTION</b>
For a U.S. specific package insert, please go to www.medacta.com or call toll free 800-901-7836.

<b>SYMBOLS / SIMBOLI / PICTOGRAMMEN / PIKTOGRAMME / PICTOGRAMAS / PICTOGRAMAS / PICTOGRAMMES / ΕΙΚΟΝΟΣΥΜΒΟΛΑ</b>	
	Do not reuse / Non riutilizzare / Niet hergebruiken / Nicht wiederverwenden / No reutilizable / Não reutilizar / Ne pas réutiliser / Μη επαναχρησιμοποιείτε
	Do not resterilize / Non ristilizzare / Niet resteriliseren / Bitte nicht re-sterilisieren / No resterilize / Não resterilizar / Ne pas restériliser / Να μην αποστειρωθεί ξανά
	Caution, read the accompanying documents / Atención, consultar i documentos d'acompanyamento / Opgelet, raadpleeg de bijgeleverde documenten / Achtung, Begleitdokumente beachten / Atención, consulte los documentos que se acompañan / Atención, consultar os documentos de acompañamento / Attention, consulter les documents d'accompagnement / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα
	Consult instructions for use / Consultare le istruzioni di utilizzo / Gebruiksaanwijzing / Bitte lesen Sie die Gebrauchsanweisung / Consulte las instrucciones antes de su uso / Consulte as instruções de uso / Consultez les instructions d'utilisation / Οδηγίες χρήσης
	Do not expose to sunlight / Non esporre alla luce del sole / Niet blootstellen aan zonlicht / Vor Sonnenlicht schützen / No exponer a la luz solar / Não expor à luz solar / Ne pas exposer à la lumière du soleil / Μη εκθέτετε σε ηλιακό φως
	Store in a dry place / Conservare in luogo asciutto / Droog bewaren / Trocken aufbewahren / Conservar en lugar seco / Conservar em local seco / Conserver au sec / Φυλάσσετε σε ξηρό χώρο
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Niet gebruiken als verpakking beschadigd is / Bitte nicht verwenden, falls Verpackung beschädigt ist / No utilizar si el envoltorio está dañado / Não utilizar se a embalagem estiver violada / Ne pas utiliser si l'emballage est abîmé / Να μην χρησιμοποιηθεί εάν η συσκευασία είναι κατεστραμμένη
	Use by / Utilizzare entro il / Te gebruiken voor / Zu verwenden bis / Utilizar antes de / A utilizar antes de / A utiliser avant / Ημερομηνία λήξης
	Lot number / Numero di lotto / Nummer van de partij / Waren-Lot-Nummer / Número de lote / Número de lote / Numéro de lot / Αριθμός παρτίδας
	Reference number / Riferimento commerciale / Handelsreferentie / Referenznummer / Numero de Referencia / Referência comercial / Référence commerciale / Νοτιέρο παραπομπής
	Sterilized by irradiation / Sterilizado mediante irradiazione / Gesteriliseerd door dostraling / Sterilisiert durch Bestrahlung / Esterilizado por irradiación / Esterilizado por irradiação / Stérilisé par irradiation / Αποστειρωμένο με ακτινοβολία
	MR Safe / Seguro per RM / MR veilig / MR geeignet / Seguro RM / Seguro para RM / RM sécurisée / ασφαλές στην μαγνητική τομογραφία

<b>ENGLISH</b>	<b>INSTRUCTIONS FOR USE MECTASCREW PEEK</b>

Important notice: the device(s) can be prescribed and implanted only by a medical doctor legally authorized to perform this type of surgery.

#### GENERAL

Before any surgery, the surgeon must be familiar with the product literature and operative technique and must read carefully these instructions for use. Patient selection is as important as implant placement or positioning. Unsuitable functional requirements may contribute to reduce the implant life. The warnings must be heeded, and the instructions for use must be strictly followed.

#### PRODUCT DESCRIPTION

These instructions for use are intended for all products described here below.

#### Medacta's product description

Medacta MectaScrew PEEK interference screws were developed to provide interference between a ligament graft and a bone tunnel in the ligament reconstruction surgery.

#### INTENDED USE

The purpose of the MectaScrew PEEK interference screws is to provide a suitable and secure intra- tunnel graft fixation in ligamentoplasty surgery.

**Indications**  
Reconstructive treatment of ruptured anterior and posterior cruciate ligaments by means of auto-, allo-, and synthetic grafts

#### Contraindications

- Osteoporosis and osteomalacia
- Adiposity or patient overweight, leading to excessive strain on the fixation
- Osteoma in the area in which the fixation button is to be placed
- Deformities of the bone, or general conditions of the bone, which preclude the implantation of the fixation button in the opinion of a physician
- Systemic diseases and metabolic disorders that may compromise the outcome of the surgery

#### WARNINGS

The success of the operation depends on compliance with the operative technique supplied, and the proper use of the dedicated instruments specially designed for that range of implants. Malpositioning may reduce implant longevity and lead to early implant failure.

One Medacta label for each device must be properly applied to the Medacta Implant Card provided, and delivered to the patient.

#### PRECAUTIONS

Ensure Medacta MectaScrew PEEK interference screws are used with the special screwdriver designed for this purpose. The screwdriver must be inserted into the screw as far as it will go in order to prevent damage to the screw and its head during insertion.

#### MRI Compatibility

The Medacta MectaScrew PEEK interference screws are MR Safe.

#### Medacta International Implants

Medacta International is not responsible for the use of its implant components in combination with a component from another manufacturer (unless otherwise specified by Medacta International in the surgical technique), therefore we advise against such a use. The components of a MectaScrew PEEK interference screw should never be reimplanted. While an implant may appear undamaged, microscopic imperfections may occur and cause implant failure. The operating surgeon must be aware that even a very small superficial damage, caused for instance by a sharp tool or electrocautering, can have an influence on the endurance of the device and can lead to fracture. Medacta MectaScrew PEEK interference screw products are non-pyrogenic.

#### Caution

The following risk factors, individually or together, may result in poor clinical outcomes:

- Inadequate bone quality, (e.g. osteoporosis, previous cruciate ligament operation)
- Systemic diseases or metabolic disorders
- History of infections or recurrent falls
- Drug dependence and abuse of alcohol and medications
- Mental incapacity of patient to understand the instructions of the physician and to comply with them
- Local bone tumors

#### INSTRUCTIONS FOR USE

#### Preoperative phase

The surgeon should verify possible patient physical limitations and mental deficiencies and he should also discuss with the patient all the details of the procedure and implant. The discussion should consider the limitations of the procedure and the constraints imposed by the selected implant. The factors which could limit the performance and stability of the implant, e.g. level of activity, should be set out to improve the patient's chances to avoid complications. The necessity to follow the postoperative instructions given by the surgeon should be fully understood by the patient. A stock of sterile implants of suitable sizes must be available and checked by the operator before surgery.

#### Handling

To avoid scratching or damaging the implants, these should be handled with the utmost care by qualified personnel and in an environment where conditions of hygiene are controlled. The implants should be kept in their undamaged packages.

#### Surgical technique

The surgeon should be fully familiar with the surgical technique. Supplementary information about the surgical techniques (brochure and video) and products are available on request. Careful preoperative planning, documented by X-rays (and/or computer tomography (CT) or magnetic resonance imaging (MRI) scans), is essential.

#### Postoperative care and follow-up

The surgeon should caution the patients to control their level of activity and avoid excessive loads on the operated limb. Moreover the surgeon should make the patients aware of the precautions to be taken in terms of exercises, treatments and limitations on activities, any limitations reported on the label, as well as exposure to magnetic fields.

Periodic follow-up and X-rays (and/or computer tomography (CT) or magnetic resonance imaging (MRI) scans) are recommended to make comparisons with the immediate postoperative condition and anticipate implant related complications. Excessive physical activity, and operated limb traumas may cause early failure of the implant. If in the case occurs, it is necessary to place the patient under supervision, evaluate the possible progression of the deterioration, and weigh the benefit of early revision.

#### ADVERSE EFFECTS AND COMPLICATIONS

Adverse effects that can occur in reconstructive treatment of ligament ruptures include:

- Infection, both deep and superficial
- Allergies, mild inflammatory and foreign body reactions to implant material
- Some adverse effects can ultimately lead to death.
- General complications include:
  - Venous thrombosis with/without pulmonary embolism
  - Cardiovascular or pulmonary disturbances
  - Haematomas
  - Systemic allergic reactions
  - Systemic pain

#### PACKAGING

The components of all MectaScrew PEEK interference screw are supplied in single-use packages.

For components delivered sterile, the sterilization method is indicated on the label. The expiration date must be checked on the label as well as the package integrity to ensure that sterility of the contents has not been compromised. If the package is damaged or has been previously opened, do not use the component. Do not resterilize.

#### STORAGE

The packages must be stored in a cool, dry place, away from light.

#### TRADEMARKS

Medacta is registered trademark of Medacta International SA, Castel San Pietro, Switzerland.

The reference text is the English text
Last Update: October 2019

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