Strada Regina  
CH - 6874 Castel San Pietro - Switzerland  
Phone +41 91 696 60 60 - Fax +41 91 696 60 66  
info@medacta.ch - www.medacta.com

Manufactured by:  
Medacta International SA  
CH-6874 Castel San Pietro - Switzerland

## MECTA TAP TI SUTURE ANCHORS

### CAUTION

For U.S. specific instructions for use, please go to [www.medacta.com](http://www.medacta.com) or call toll free 800-901-7836.

#### SYMBOLS / SIMBOLI / PICTOGRAMMES / PIKTOGRAMME / PICTOGRAMAS / PICTOGRAMAS / PICTOGRAMMES / ΕΙΚΟΝΟΣΥΜΒΟΛΑ

	Do not reuse / Non riutilizzare / Niet hergebruiken / Nicht wiederverwenden / No reutilizable / Não reutilizar / Ne pas réutiliser / Μην επαναχρησιμοποιείτε
	Do not resterilize / Non risterilizzare / Niet resteriliseren / Bitte nicht re-sterilisieren / No resterilize / Não resterilizar / Ne pas restériliser / Να μην αποστειρωθεί ξανά
	Caution, read the accompanying documents / Attenzione, consultare i documenti di accompagnamento / Opgelet, raadpleeg de bijgeleverde documenten / Achtung, Begleitdokumente beachten / Atención, consulte los documentos que se acompañan / Atenção, consultar os documentos de acompanhamento / Attention, consulter les documents d'accompagnement / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα
	Consult instructions for use / Consultare le istruzioni di utilizzo / Gebruiksaanwijzing / Bitte lesen Sie die Gebrauchsanweisung / Consulte las instrucciones antes de su uso / Consulte as instruções de uso / Consultez les instructions d'utilisation / Οδηγίες χρήσης
	Do not expose to sunlight / Non esporre alla luce del sole / Niet blootstellen aan zonlicht / Vor Sonnenlicht schützen / No exponer a la luz solar / Não expor à luz solar / Ne pas exposer à la lumière du soleil / Μην εκθέτετε σε ηλιακό φως
	Store in a dry place / Conservare in luogo asciutto / Droog bewaren / Trocken aufbewahren / Conservar en lugar seco / Conservar em local seco / Conserver au sec / Φυλάσσετε σε ξηρό χώρο
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Niet gebruiken als verpakking beschadigd is / Bitte nicht verwenden, falls Verpackung beschädigt ist / No utilize si el envoltorio está dañado / Não utilizar se a embalagem estiver violada / Ne pas utiliser si l'emballage est abîmé / Να μην χρησιμοποιηθεί εάν η συσκευασία είναι κατεστραμμένη
	Use by / Utilizzare entro il / Te gebruiken voor / Zu verwenden bis / Utilizar antes de / A utilizar antes de / A utiliser avant / Ημερομηνία λήξης
	Lot number / Numero di lotto / Nummer van de partij / Waren-Lot-Nummer / Número de lote / Número de lote / Numéro de lot / Αριθμός παρτίδας
	Reference number / Riferimento commerciale / Handelsreferentie / Referenznummer / Numero de Referencia / Referência comercial / Référence commerciale / Νούμερο παραπομπής
	Sterilized with ethylene oxide / Sterilizzato con ossido di etilene / Gesteriliseerd met ethyleenoxide / Ethylenoxid-sterilisiert / Esterilizado con óxido de etileno / Esterilizado com óxido de etileno / Stérilisé à l'oxyde d'éthylène / Αποστειρωμένο με οξείδιο του αιθυλενίου
	MR Conditional / A compatibilità condizionata con la risonanza magnetica / MR-conditioneel / Bedingt MR tauglich / Compatibles con RM bajo ciertas condiciones / Qualificado para ressonância magnética / RM conditionnelle / Μαγνητικός συντονισμός υπό προϋποθέσεις

## Index

ENGLISH: MECTATAP TI SUTURE ANCHOR - INSTRUCTIONS FOR USE .....	3
ITALIANO: MECTATAP TI SUTURE ANCHOR - ISTRUZIONI PER L'USO .....	8
NEDERLANDS: MECTATAP TI SUTURE ANCHOR - GEBRUIKSINSTRUCTIES.....	13
DEUTSCH: MECTATAP TI SUTURE ANCHOR - GEBRAUCHSANWEISUNGEN .....	18
ESPAÑOL: MECTA AP TI SUTURE ANCHOR - INSTRUCCIONES DE USO .....	23
PORTUGUÊS: MECTATAP TI SUTURE ANCHOR - INSTRUÇÕES DE UTILIZAÇÃO .....	28
FRANÇAIS: MECTATAP TI SUTURE ANCHOR - INSTRUCTIONS D'UTILISATION .....	33
ΕΛΛΗΝΙΚΑ: MECTATAP TI SUTURE ANCHOR - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ.....	38

## ENGLISH: MECTATAP TI SUTURE ANCHOR - INSTRUCTIONS FOR USE

Important notice: the device(s) can be prescribed and implanted only by a medical doctor legally authorized to perform this type of surgery.

### 1. GENERAL

Before any surgery, the surgeon must be familiar with the product literature and operative technique and must read carefully these instructions for use. Patient selection is as important as implant placement or positioning. Unsuitable functional requirements may contribute to reduce the implant life. The warnings must be heeded, and the instructions for use must be strictly followed.

### 2. PRODUCT DESCRIPTION

These instructions for use are intended for all products described here below.

CATALOGUE NUMBER	ANCHOR DESCRIPTION
05.10.013	MectaTap TI Ø5.0
05.10.014	MectaTap TI Ø6.5

Medacta MectaTap TI Suture Anchors are manufactured in Titanium alloy.

Medacta MectaTap TI Suture Anchors are provided with sutures directly preloaded into the implant.

CATALOGUE NUMBER	INCLUDED SUTURES
05.10.013	
05.10.014	2x USP#2/EP#5 HS Fiber 39» preloaded sutures

### Material Specifications

Anchors: Ti6Al4V ELI ISO 5832-3 [MectaTap TI]

Driver shaft/body: Stainless steel

Driver Handle: Polycarbonate medical grade

Suture: non-absorbable USP#2/EP#5 HS Fiber 39» Ultra-High molecular weight polyethylene (average diameter: 0.500-0.599mm)

### 3. INTENDED USE / INDICATIONS

The MectaTap TI Suture anchors are intended to be used for soft tissue refixation within the shoulder joint (e.g.: cuff rotator repair, biceps tenodesis).

The device may be used in either arthroscopic or open surgical procedures. After the suture is anchored to the bone by the anchor, it may be used to reattach soft tissue to the bone. The device, in that way, stabilizes the damaged soft tissues, in combination with appropriate postoperative immobilization, throughout the healing period.

#### 4. CONTRAINDICATIONS

- Pathological conditions in the soft tissue to be repaired or reconstructed which would adversely affect suture fixation.
- Pathological conditions of bone which adversely affect stability of the MectaTap TI Suture Anchors.
- Attachment of artificial ligament or other implants.
- Foreign body sensitivity, known or suspected allergies to implant and/or instrument materials.

#### 5. WARNINGS AND PRECAUTIONS

The success of the operation depends on compliance with the operative technique supplied, and the proper use of the dedicated instruments especially designed for that range of implants. Surgeon must choose proper implant size based on specific procedure and patient history. Malpositioning may reduce implant longevity and lead to early implant failure. The disposable driver is single use and must be disposed of according to hospital policy and procedure.

One Medacta label for each device must be properly applied to the Medacta Implant Card provided, and delivered to the patient.

The risk of MectaTap TI Suture Anchor breakage during implantation is reduced by following the specified instructions for use listed below. Proper selection and placement of the implant are important considerations in the successful utilization of this device. Proper orientation and alignment of instruments is important during implantation of the MectaTap TI Suture Anchor to minimize the risk of breakage of the anchor. Breakage of the MectaTap TI Suture Anchor can occur if:

- Any bending or shear load is applied on the anchor/insertor while inserting the MectaTap TI anchor;
- The MectaTap TI Suture Anchor driver is used for prying.

Do not disengage the sutures from the disposable driver while inserting the MectaTap TI Suture Anchor before the final positioning depth is reached. Do not use any other suture other than non-absorbable USP#2/EP#5 UHMWPE suture.

#### MRI Compatibility

Non-clinical testing has demonstrated titanium implants are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions. Static magnetic field of 1.5 Tesla and 3 Tesla, with :

- Maximum spatial field gradient of 12,800 G/cm (128 T/m)
- Maximum force product of 231,000,000 G<sup>2</sup>/cm (231 T<sup>2</sup>/m)

Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode). Under the scan conditions defined above, titanium implants expected to produce a maximum temperature rise of less than:

- 1.5°C (2 W/kg, 1.5 Tesla) RF-related temperature increase with a background temperature increase of ≈1.3°C (2 W/kg, 1.5 Tesla)

- 3.0°C (2 W/kg, 3 Tesla) RF-related temperature increase with a background temperature increase of ≈2.3°C (2 W/kg, 3 Tesla)

after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 20.5 mm from the titanium implant when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

## 6. MEDACTA INTERNATIONAL IMPLANTS

Medacta International is not responsible for the use of its implant components in combination with a component from another manufacturer (unless otherwise specified by Medacta International in the surgical steps), therefore we advise against such practice. The MectaTap TI Suture Anchors should never be re-implanted. While an implant may appear undamaged, microscopic imperfections may occur and cause implant failure. The operating surgeon must be aware that even a very small superficial damage, caused for instance by a sharp tool or electrocautery, can have an influence on the endurance of the device and can lead to failure.

MectaTap TI Suture Anchors are non-pyrogenic.

## 7. RISK FACTORS

The following risk factors, individually or together, may result in poor clinical outcomes:

- Inadequate bone quality, (e.g. osteoporosis, osteopenia)
- Systemic diseases or metabolic disorders
- History of infections or recurrent falls
- Drug dependence and abuse of alcohol and medicaments
- Mental incapacity of patient to understand the instructions of the physician and to comply with them
- Local bone tumors
- Physical condition that would eliminate, or tend to eliminate, adequate implant support or impede healing
- Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period

## 8. INSTRUCTIONS FOR USE

### PREOPERATIVE PHASE

The surgeon should verify possible patient physical limitations and mental deficiencies. The surgeon should also discuss with the patient all the details of the procedure and implant. The discussion should consider the limitations of the procedure and the constraints imposed by the selected implant. The factors which could limit the performance and stability of the implant, (e.g. level of activity), should be set out to improve the patient's chances to avoid complications. The necessity to follow the postoperative instructions given by the surgeon should be fully understood by the patient. A stock of sterile implants of suitable sizes must be available and checked by the operator before surgery.

## 9. HANDLING


To avoid scratching or damaging the implants, these must be handled with the utmost care by qualified personnel. Sterile implants must be opened with the use of sterile technique in a controlled environment. The implants must be kept in their undamaged packages.

## 10. SURGICAL TECHNIQUE


The surgeon should be fully familiar with the surgical technique. Supplementary information about the surgical techniques (brochure and video) and products are available on request. Careful preoperative planning is essential.

### MectaTap TI Suture Anchors:

- Find the desired location for the placement of the chosen MectaTap TI Suture Anchor; insert the preloaded driver with the MectaTap TI Suture Anchor through the proper arthroscopic portal to reach the desired location
- Gently tap on the back handle of the disposable driver to let the MectaTap TI Suture Anchor's tip pierce the bone in the selected spot
- Keep tapping with light and repeated strokes on the back handle of the disposable driver to let the MectaTap TI Suture Anchor enter completely into the bone until the distal laser marking on the driver is flush with the bone surface

- Unlock the sutures from the plastic handle and gently disengage the driver's tip from the MectaTap TI Suture Anchor, pulling the driver out through the arthroscopic portal

- Manage the 2x USP#2/EP#5 HS Fiber 39" UHMWPE sutures through soft tissue with surgeon's preferred suture-passing and knotting technique
- Palpate the repaired soft tissue using an arthroscopic probe to ensure that the repair is secure
- Cut and remove suture in excess

### 11. POSTOPERATIVE CARE AND FOLLOW-UP

The surgeon should caution the patients to control their level of activity and avoid excessive loads on the operated limb. Moreover the surgeon should make the patients aware of the precautions to be taken in terms of exercises, treatments and limitations on activities, any limitations reported on the label, as well as exposure to magnetic fields. The patient must be told that implants can affect the results of computer tomography (CT) or magnetic resonance imaging (MRI) scans. Periodic follow-up visits are recommended to make comparisons with the immediate postoperative condition and anticipate implant related complications. Excessive physical activity, and operated limb traumas may cause early failure of the implant. If the case occurs, it is necessary to place the patient under supervision, evaluate the possible progression of the deterioration, and weigh the benefit of early revision.

## 12. ADVERSE EFFECTS AND COMPLICATIONS

Adverse effects that can occur in the reattachment of soft tissue to bone in orthopedic surgical procedures:

- Infection, both deep and superficial
- Allergies, mild inflammatory and foreign body reactions to implant material

Some adverse effects can ultimately lead to death.

General complications include:

- Venous thrombosis with/without pulmonary embolism;
- Cardiovascular or pulmonary disturbances;
- Haematomas;
- Systemic allergic reactions;
- Systemic pain.

## 13. PACKAGING

The Mectatap TI Suture Anchors are supplied in single-use packages. For components delivered sterile, the sterilization method is indicated on the label. The expiration date must be checked on the label as well as the package integrity to ensure that sterility of the contents has not been compromised. If the package is damaged or has been previously opened, do not use the component. Do not re-sterilize.

## 14. STORAGE

The packages must be stored in a cool, dry place, away from light.

## 15. TRADEMARKS

Medacta® is registered trademark of Medacta International SA, Castel San Pietro, Switzerland.

The reference text is the English text.

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