Strada Regina
CH - 6874 Castel San Pietro - Switzerland
Phone +41 91 696 60 60 - Fax +41 91 696 60 66
info@medacta.ch - www.medacta.com


Manufactured by:
Medacta International SA
CH-6874 Castel San Pietro - Switzerland

MEDACTA SUTURE FAMILY - POWERKNOT HIGH STRENGTH SUTURE

CAUTION

For U.S. specific instructions for use, please go to www.medacta.com or call toll free 800-901-7836.

SYMBOLS / SIMBOLI / PICTOGRAMMES / PIKTOGRAMME / PICTOGRAMAS / PICTOGRAMAS / PICTOGRAMMES / ΕΙΚΟΝΟΣΥΜΒΟΛΑ

	Do not reuse / Non riutilizzare / Niet hergebruiken / Nicht wiederverwenden / No reutilizable / Não reutilizar / Ne pas réutiliser / Μην επαναχρησιμοποιείτε
	Do not resterilize / Non risterilizzare / Niet resteriliseren / Bitte nicht re-sterilisieren / No resterilize / Não resterilizar / Ne pas restériliser / Να μην αποστειρωθεί ξανά
	Caution, read the accompanying documents / Attenzione, consultare i documenti di accompagnamento / Opgelet, raadpleeg de bijgeleverde documenten / Achtung, Begleitdokumente beachten / Atención, consulte los documentos que se acompañan / Atenção, consultar os documentos de acompanhamento / Attention, consulter les documents d'accompagnement / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα
	Consult instructions for use / Consultare le istruzioni di utilizzo / Gebruiksaanwijzing / Bitte lesen Sie die Gebrauchsanweisung / Consulte las Instrucciones antes de su uso / Consulte as instruções de uso / Consultez les instructions d'utilisation / Οδηγίες χρήσης
	Do not expose to sunlight / Non esporre alla luce del sole / Niet blootstellen aan zonlicht / Vor Sonnenlicht schützen / No exponer a la luz solar / Não expor à luz solar / Ne pas exposer à la lumière du soleil / Μην εκθέτετε σε ηλιακό φως
	Store in a dry place / Conservare in luogo asciutto / Droog bewaren / Trocken aufbewahren / Conservar en lugar seco / Conservar em local seco / Conserver au sec / Φυλάσσετε σε ξηρό χώρο
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Niet gebruiken als verpakking beschadigd is / Bitte nicht verwenden, falls Verpackung beschädigt ist / No utilize si el envoltorio está dañado / Não utilizar se a embalagem estiver violada / Ne pas utiliser si l'emballage est abîmé / Να μην χρησιμοποιηθεί εάν η συσκευασία είναι κατεστραμμένη
	Use by / Utilizzare entro il / Te gebruiken voor / Zu verwenden bis / Utilizar antes de / A utilizar antes de / A utiliser avant / Ημερομηνία λήξης
	Lot number / Numero di lotto / Nummer van de partij / Waren-Lot-Nummer / Número de lote / Número de lote / Numéro de lot / Αριθμός παρτίδας
	Reference number / Riferimento commerciale / Handelsreferentie / Referenznummer / Numero de Referencia / Referência comercial / Référence commerciale / Νούμερο παραπομπής
	Sterilized with ethylene oxide / Sterilizzato con ossido di etilene / Gesteriliseerd met ethyleenoxide / Ethylenoxid-sterilisiert / Esterilizado con óxido de etileno / Esterilizado com óxido de etileno / Stérilisé à l'oxyde d'éthylène / Αποστειρωμένο με οξειδίο του αιθυλενίου
	Quantity / Quantità / Hoeveelheid / Quantität / Cantidad / Quantidade / Quantité / Ποσότητα
	Manufacturer / Produttore / Fabrikant / Hersteller / Fabricante / Fabricante / Fabricant / Κατασκευαστής
	MR Safe / Compatible con risonanza magnetica / MR-veilig / MR-sicher / Seguro ante resonancia magnética / Seguro para ressonância magnética / Seguro para ressonância magnética / Ασφαλές στην μαγνητική τομογραφία

Index

ENGLISH: POWERKNOT - INSTRUCTIONS FOR USE	3
ITALIANO: POWERKNOT - ISTRUZIONI PER L'USO	5
NEDERLANDS: POWERKNOT - GEBRUIKSINSTRUCTIES.....	7
DEUTSCH: POWERKNOT - GEBRAUCHSANWEISUNGEN	9
ESPAÑOL: POWERKNOT - INSTRUCCIONES DE USO	11
PORTUGUÊS: POWERKNOT - INSTRUÇÕES DE UTILIZAÇÃO	13
FRANÇAIS: POWERKNOT - INSTRUCTIONS D'UTILISATION.....	15
ΕΛΛΗΝΙΚΑ: POWERKNOT - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	17

ENGLISH: POWERKNOT - INSTRUCTIONS FOR USE

Important notice: the device(s) should only be prescribed and implanted by a medical doctor legally authorized to perform this type of surgery.

1. GENERAL

Before any surgery, the surgeon must be familiar with the product literature and operative technique and must read carefully these instructions for use. Patient selection is as important as implant placement or positioning. Unsuitable functional requirements may contribute to reduce the implant life. The warnings must be heeded and the instructions for use must be strictly followed.

2. PRODUCT DESCRIPTION

Medacta Suture Family is comprised of non-absorbable sutures made of 100% Ultra High Molecular Weight Polyethylene (UHMWPE). PowerKnot and PowerKnot RDI High Strength Suture is a braided suture made of UHMWPE over a core of UHMWPE yarns. PowerKnot RDI High Strength Suture comprises a color asymmetric pattern providing a directional indication on the outer surface of the suture, aiding the surgeon when applying the suture. The sutures meet or exceed U.S. Pharmacopeia standards for non-absorbable surgical sutures (except for diameter requirements). The diameters of the sutures typically span the diameter requirements of two USP sizes. For that reason, Medacta tests the suture to conform to the greater USP size specification requirement. The Medacta Suture Family is supplied sterile in a range of gauge sizes and pre-cut lengths with heat stiffened ends, and in some cases, with various curved needles made of stainless steel. The Medacta Suture family comprises non-dyed and dyed yarns. Dyes may include: chromium-cobalt-aluminum oxide spinel for the blue colour, high purity carbon for the black colour and bismuth trioxide (Bi₂O₃) for the yellow colour.

3. INTENDED USE / INDICATIONS

Medacta Suture Family is indicated for soft tissue approximation and/or ligation. These sutures may be incorporated, as components, into surgeries where constructs including those with allograft or autograft tissues are used for repair.

4. CONTRAINDICATIONS

Medacta Suture Family is not for use in cardiac indications

5. WARNINGS AND PRECAUTIONS

Check the expiration date prior to use. Medacta Family Sutures are meant to be single use. Reutilizing them is strictly forbidden and would lead to unpredictable outcomes. Do not touch sharp areas with fingers. Injury may result. Visually inspect the devices prior and after use in order to verify that they did not experience any mechanical damage. Any non-functional device should be immediately returned to Medacta® along with its label. The type of malfunction should be reported.

MRI Compatibility

Medacta Suture Family has been evaluated for safety and compatibility in the MR environment, and this device is not likely to result in patient harm during MR scanning.

6. MEDACTA INTERNATIONAL IMPLANTS

Medacta Suture Family should never be re-implanted. While an implant may appear undamaged, microscopic imperfections may occur and cause implant failure. The operating surgeon must be aware that even a very small superficial damage, caused for instance by a sharp tool, a radiofrequency device or electrocautery, can have an influence on the endurance of the device and can lead to failure.

7. PRECAUTIONS

In handling Medacta Suture products, care should be taken to avoid damage from handling. Pay attention to crushing and crimping derived from the application of forceps or needle holders. Avoid grasping the needle at the point or swage to avoid damage of these areas. Reshaped needles may lose penetration performance, strength and resistance to bending and breaking. Dispose according to hospital guidelines pertaining to sharp disposable instruments. Make sure that all the knots have been secured using recognized surgical knot tying techniques.

8. INSTRUCTIONS FOR USE

Knot tying requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstance and the experience of the final user. Care should be exercised during handling and use to avoid suture damage. Precautions should be taken when using surgical instruments to grasp the suture. Instruments used should be in good working condition with no sharp edges to avoid crushing or severing of the suture. The suture should only be held at the free end.

9. ADVERSE EFFECTS AND COMPLICATIONS

Common non-absorbable suture reactions may include:

- Wound dehiscence and erythema at the wound site
- Inflammatory tissue reaction
- Enhanced bacterial infectivity
- Calculi formation in urinary and biliary tracts
- Edema
- Pain

10. PACKAGING

The Medacta Suture Family is provided disposable and sterile and is supplied in boxes of multiple units each. The sterilization method used is indicated on the label. The expiration date must be checked on the label as well as the package integrity to ensure that sterility of the contents has not been compromised. If the package is damaged or has been previously opened, do not use the components.

11. STORAGE AND HANDLING

The packages must be stored in a cool, dry place, away from light. Handle with care.

12. TRADEMARKS

Medacta® is a registered trademark of Medacta International SA, Castel San Pietro, Switzerland.

The reference text is the English text.

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