Bo.Ce.Mix femoral blood sponge

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Bo.Ce.Mix femoral blood sponge

Read these instructions carefully before use

This Medical Device was developed and manufactured under quality systems certified as conforming to the EN ISO 9001 – EN ISO 13485 standards and to conform with the essential requirements of Directive 93/42/ECC, 2007/47/EC (and subsequent amendments to these directives).

The device is classified as Class IIa sterile according to annex IX, rule 6 of the directive. It is intended for SINGLE-USE or SINGLE-PATIENT USE (for repeated short term transient use on the same patient) and is sterilised with Ethylene Oxide gas.

Improper use may lead to a risk of cross-infection or spreading of disease, with serious repercussions to patient health.

We are not responsible for the improper or incorrect use of the Device, or for lack of control before use.

Reusing the device is strictly prohibited, as it can seriously affect patient's health, since it is not possible to accurately monitor the possibility of thoroughly cleaning the device or to ensure that the device has been effectively re-sterilized.

Before using:

- 5) Please ensure that device is within its expiry date (within 5 years of the date of sterilization)
- 6) Please ensure that the packaging and carton are in good condition, are not wet, and have not been opened, tampered with or altered.
- 7) Check that the device is not broken or faulty and that all parts appear clean and in good working order.
- 8) If any of the previous conditions apply, **do not** use the Device.

Instructions for Use:

- 4) Open the packet containing the device, remove it for inspection and, if necessary, lay it on a towel or sterile surface.
- 5) Connect the connection funnel to the suction system and, through the control hole, adjust the vacuum according to the specific needs.

After use, the Device can be eliminated/disposed of in accordance with local hospital waste regulations. The plastic used to make Bo.Ce.Mix femoral blood sponge is non-toxic and can be incinerated without producing toxic gasses or other harmful or polluting elements.

= Contains phthalates [Bis(2-ethylhexyl)phthalate]

Check that the labels are intact, are not damp, and have not deteriorated, that the data on the packagings and the instructions for use have not been tampered with in any way.

If any of the above mentioned conditions apply, please contact the Manufacturer before use.