MECTALOCK PEEK SUTURE ANCHORS MECTALOCK C SUTURE ANCHORS MECTALOCK TI SUTURE ANCHORS

CAUTION

For a U.S. specific package insert, please go to www.medacta.com or call toll free 800-901-7836.

SYMBOLS / SIMBOLI / PICTOGRAMMES / PIKTOGRAMME / PICTOGRAMAS / PICTOGRAMAS / PICTOGRAMMES / ΕΙΚΟΝΟΣΥΜΒΟΛΑ

APPLIES TO ALL DEVICES

	Do not reuse / Non riutilizzare / Niet hergebruiken / Nicht wiederverwenden / No reutilizable / Não reutilizar / Ne pas réutiliser / Μην επαναχρησιμοποιείτε
	Do not resterilize / Non risterilizzare / Niet resteriliseren / Bitte nicht re-sterilisieren / No reesterilize / Não reesterilizar / Ne pas résteriliser / Na μην αποστειρωθεί ξανά
	Caution, read the accompanying documents / Attenzione, consultare i documenti di accompagnamento / Opgelet, raadpleeg de bijgeleverde documenten / Achtung, Begleitdokumente beachten / Atención, consulte los documentos que se acompañan / Atenção, consultar os documentos de acompanhamento / Attention, consulter les documents d'accompagnement / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα
	Consult instructions for use / Consultare le istruzioni di utilizzo / Gebruiksaanwijzing / Bitte lesen Sie die Gebrauchsanweisung / Consulte las Instrucciones antes de su uso / Consulte as instruções de uso / Consultez les instructions d'utilisation / Οδηγίες χρήσης
	Do not expose to sunlight / Non esporre alla luce del sole / Niet blootstellen aan zonlicht / Vor Sonnenlicht schützen / No exponer a la luz solar / Não expor à luz solar / Ne pas exposer à la lumière du soleil / Μην εκθέτετε σε ηλιακό φως
	Store in a dry place / Conservare in luogo asciutto / Droog bewaren / Trocken aufbewahren / Conservar en lugar seco / Conservar em local seco / Conserver au sec / Φυλάσσετε σε ξηρό χώρο
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Niet gebruiken als verpakking beschadigd is / Bitte nicht verwenden, falls Verpackung beschädigt ist / No utilice si el envoltorio está dañado / Não utilizar se a embalagem estiver violada / Ne pas utiliser si l'emballage est abîmé / Na μην χρησιμοποιηθεί εάν η συσκευασία είναι κατεστραμμένη
	Use by / Utilizzare entro il / Te gebruiken voor / Zu verwenden bis / Utilizar antes de / A utilizar antes de / A utiliser avant / Ημερομηνία λήξης
	Lot number / Numero di lotto / Nummer van de partij / Waren-Lot-Nummer / Número de lote / Número de lote / Numéro de lot / Αριθμός παρτίδας
	Reference number / Riferimento commerciale / Handelsreferentie / Referenznummer / Numero de Referencia / Referência comercial / Référence commerciale / Νούμερο παραπομπής
	Sterilized with ethylene oxide / Sterilizzato con ossido di etilene / Gesteriliseerd met ethyleenoxide / Ethylenoxid-sterilisiert / Esterilizado con óxido de etileno / Esterilizado com óxido de etileno / Stérilisé à l'oxyde d'éthylène / Αποστειρωμένο με οξείδιο του αιθυλενίου
	Sterilized by irradiation / Sterilizzato mediante irradiazione / Gesteriliseerd door doorstraling / Sterilisiert durch Bestrahlung / Esterilizado por irradiación / Esterilizado por irradiação / Stérilisé par irradiation / Αποστειρωμένο με ακτινοβολία
	Manufacturer / Produttore / Fabrikant / Hersteller / Fabricante / Fabricante / Fabricant / Κατασκευαστής
	MR Conditional / A compatibilità condizionata con la risonanza magnetica / MR-conditioneel / Bedingt MR tauglich / Compatibles con RM bajo ciertas condiciones / Qualificado para ressonância magnética / RM conditionnelle / Μαγνητικός συντονισμός υπό προϋποθέσεις
	MR Safe / Sicuro per RM / MR veilig / MR geeignet / Seguro para RM / Seguro para RM / Compatible avec IRM / Ασφαλές στην μαγνητική τομογραφία
	Temperature limit
	Single sterile barrier system / Sistema a Singola Barriera Sterile / Enkele steriele verpakking / Einfaches steriles Barriersystem / Sistema de barrera estéril única / Sistema de Barreira Estéril Única / Système à barrière stérile simple / Σύστημα μονού στείρου εμποδίου
	Double sterile barrier system / Sistema a Doppia Barriera Sterile / Dubbele steriele verpakking / Doppeltes steriles Barriersystem / Sistema de doble barrera estéril / Sistema de Dupla Barreira Estéril / Système à barrière stérile double / Σύστημα διπλού στείρου εμποδίου

Index

ENGLISH: MECTALOCK SUTURE ANCHORS - INSTRUCTIONS FOR USE	3
ITALIANO: ANCORE PER SUTURA MECTALOCK - ISTRUZIONI PER L'USO.....	11
NEDERLANDS: HANDLEIDING MECTALOCK HECHTANKERS - GEBRUIKSINSTRUCTIES.....	19
DEUTSCH: MECTALOCK FADENANKER - GEBRAUCHSANWEISUNG ZUR VERWENDUNG	27
ESPAÑOL: MECTALOCK SUTURE ANCHORS - INSTRUCCIONES DE USO	35
PORTUGUÊS: MECTALOCK SUTURE ANCHORS - INSTRUÇÕES DE UTILIZAÇÃO	44
FRANÇAIS: MECTALOCK SUTURE ANCHORS - INSTRUCTIONS D'UTILISATION.....	52
ΕΛΛΗΝΙΚΑ: MECTALOCK SUTURE ANCHORS - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	60

ENGLISH: MECTALOCK SUTURE ANCHORS - INSTRUCTIONS FOR USE

Important notice: the device(s) can be prescribed and implanted only by a doctor legally authorized to perform this type of surgery.

1. GENERAL

Before any surgery, the surgeon must be familiar with the sales product literature and operative technique and must carefully read these instructions for use. Patient selection is as important as implant placement or positioning. Overweight patients or unsuitable functional requirements may generate exceptional stresses and reduce the implant life. The warnings must be heeded, and the instructions for use must be strictly followed

2. PRODUCT DESCRIPTION

These instructions for use are intended for all products described here below.

CATALOGUE NUMBER	ANCHOR DESCRIPTION	COMPATIBLE REUSABLE DRILL/PUNCH
05.11.001	MectaLock PEEK Ø2.4 Long driver	05.11.10.0035
05.11.101	MectaLock C Ø2.4 Long driver	
05.11.002	MectaLock PEEK Ø2.9 Long driver	05.11.10.0036
05.11.102	MectaLock C Ø2.9 Long driver	
05.10.001	MectaLock PEEK Ø2.4 Short driver	05.10.10.0015
05.10.101	MectaLock C Ø2.4 Short driver	
05.10.002	MectaLock PEEK Ø2.9 Short driver	05.10.10.0016
05.10.102	MectaLock C Ø2.9 Short driver	
05.10.004	MectaLock PEEK Ø4.5 Short driver	05.10.10.0021
05.10.104	MectaLock C Ø4.5 Short driver	
05.10.005	MectaLock PEEK Ø5.5 Short driver	05.10.10.0022
05.10.105	MectaLock C Ø5.5 Short driver	
05.10.006	MectaLock PEEK Ø6.5 Short driver	05.10.10.0023
05.10.106	MectaLock C Ø6.5 Short driver	
05.10.011	MectaLock TI Ø5.0	Not needed
05.10.012	MectaLock TI Ø6.5	Not needed

Medacta MectaLock PEEK and MectaLock C (knotless) and MectaLock TI (preloaded) are manufactured respectively in PolyEtherEtherKeton (PEEK), biocomposite or Titanium alloy. The MectaLock Suture Anchors are provided with sutures directly preloaded into the implant (MectaLock TI) or placed directly into the final packaging (MectaLock PEEK and MectaLock C).

CATALOGUE NUMBER	ANCHOR FAMILY	INCLUDED SUTURES
05.10.004	MectaLock PEEK	1x USP#2/EP#5 HS Fiber 39" UHMWPE suture
05.10.005		
05.10.006		
05.10.001		
05.10.002		
05.11.001		
05.11.002		
05.10.101	MectaLock C	1x USP#2/EP#5 HS Fiber 39" UHMWPE suture
05.10.102		
05.10.104		
05.10.105		
05.10.106		
05.11.101		
05.11.102		
05.10.011	MectaLock TI	2x USP#2/EP#5 HS Fiber 39" preloaded sutures
05.10.012		

Anchors: PolyEtherEtherKeton (PEEK) [MectaLock PEEK]; Ti6Al4V ELI ISO 5832-3 [MectaLock TI] and RESOMER® LR 706 S + β -TCP (MectaLock C)

Driver shaft/body: Stainless steel

Driver Handle: Polycarbonate medical grade

Suture: non-absorbable USP#2/EP#5 HS Fiber 39" Ultra-High molecular weight polyethylene (average diameter: 0.500-0.599mm (MectaLock PEEK & MectaLock TI) and non-absorbable USP#2/EP#5 HS Fiber 39" UHMWPE suture (MectaLock C)

3. INTENDED USE / INDICATIONS

The device is indicated for use in both arthroscopic and open surgical procedures. Soft tissues hereby reattached to bone, including the following system specific indications:

- MectaLock TI Suture anchors: soft tissue refixation within the shoulder joint (i.e.: rotator cuff repair, biceps tenodesis)
- MectaLock PEEK & MectaLock C Suture anchors size 2.4 & 2.9 mm: treatment of hip and shoulder instability by means of refixation of the acetabular or glenoid labrum to the bone.
- MectaLock PEEK & MectaLock C Suture anchors size 4.5, 5.5 & 6.5 mm: soft tissue refixation within the shoulder joint (i.e.: rotator cuff repair)

The suture anchor system in that way stabilizes the damaged soft tissue, in combination with appropriate postoperative immobilization, throughout the healing period.

4. CONTRAINDICATIONS

- Pathological conditions in the soft tissue to be repaired or reconstructed which would adversely affect suture fixation
- Pathological conditions of bone which adversely affect the MectaLock Suture Anchors
- Physical condition that would eliminate, or tend to eliminate, adequate implant support or retard healing
- Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period
- Attachment of artificial ligament or other implants
- Foreign body sensitivity, known or suspected allergies to implant and/or instrument materials

5. WARNINGS AND PRECAUTIONS

The success of the operation depends on compliance with the operative technique supplied, and the proper use of the dedicated instruments specially designed for that range of implants. Surgeon must choose proper implant size based on specific procedure and patient history. Malpositioning may reduce implant longevity and lead to early implant failure. The disposable driver is single use and must be disposed of according to hospital policy and procedure.

One Medacta label for each device must be properly applied to the Medacta Implant Card provided and delivered to the patient.

Ensure the MectaLock Suture Anchors are used with the corresponding compatible reusable drill bit or punch. Use the corresponding compatible reusable drill bit to create the pilot hole or a punch to create the pilot socket. Maintain proper alignment during insertion of the anchor and prevent disengagement from the driver. The risk of MectaLock Suture Anchor breakage during implantation is reduced by following the specified instructions for use listed below. Proper selection and placement of the implant are important considerations in the successful utilization of this device. Proper orientation and alignment of instruments is important during implantation of the MectaLock Suture Anchor to minimize possible breakage of the anchor. Breakage of the MectaLock Suture Anchor can occur if:

- The pilot hole is not drilled to the proper depth
- The MectaLock Suture Anchor is not properly aligned with the pilot hole
- Any lateral loading is applied on the anchor/insertor while inserting the MectaLock anchor
- The MectaLock Suture Anchor driver is used as a lever

Do not disengage the sutures from the disposable driver while inserting the MectaLock TI Suture Anchor before the final positioning depth is reached. Do not use any other suture other than non-absorbable USP#2/EP#5 UHMWPE suture.

MRI Compatibility

Non-clinical testing has demonstrated titanium implants are MR Conditional (MectaLock TI). A patient with this device can be safely scanned in an MR system meeting the following conditions. Static magnetic field of 1.5 Tesla and 3 Tesla, with:

- Maximum spatial field gradient of 12,800 G/cm (128 T/m)
- Maximum force product of 231,000,000 G²/cm (231 T²/m)

Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode). Under the scan conditions defined above, titanium implants expected to produce a maximum temperature rise of less than:

- 1.5°C (2 W/kg, 1.5 Tesla) RF-related temperature increase with a background temperature increase of ≈1.3°C (2 W/kg, 1.5 Tesla)
- 3.0°C (2 W/kg, 3 Tesla) RF-related temperature increase with a background temperature increase of ≈2.3°C (2 W/kg, 3 Tesla)

after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 20.5 mm from the titanium implant when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

The MectaLock PEEK and MectaLock C Suture Anchors are MR Safe.

6. MEDACTA INTERNATIONAL IMPLANTS

Medacta International is not responsible for the use of its implant components in combination with a component from another manufacturer (unless otherwise specified by Medacta International in the surgical steps), therefore we advise against such a use. The MectaLock Suture Anchors should never be reimplanted. While an implant may appear undamaged, microscopic imperfections may occur and cause implant failure. The operating surgeon must be aware that even a very small superficial damage, caused for instance by a sharp tool or electrocautery, can have an influence on the endurance of the device and can lead to failure. Medacta MectaLock Suture Anchors are non-pyrogenic.

7. RISK FACTORS

The following risk factors, individually or together, may result in poor clinical outcomes:

- Inadequate bone quality, (e.g. osteoporosis, osteopenia)
- Systemic diseases or metabolic disorders
- History of infections or recurrent falls
- Drug dependence and abuse of alcohol and medicaments
- Mental incapacity of patient to understand the instructions of the physician and to comply with them
- Local bone tumors

8. INSTRUCTIONS FOR USE

PREOPERATIVE PHASE

The surgeon should verify possible patient physical limitations and mental deficiencies and he should also discuss with the patient all the details of the procedure and implant. The discussion should consider the limitations of the procedure and the constraints imposed by the selected implant. The factors which could limit the performance and stability of the implant, e.g. level of activity, should be set out to improve the patient's chances to avoid complications. The necessity to follow the postoperative instructions given by the surgeon should be fully understood by the patient. A stock of sterile implants of suitable sizes must be available and checked by the operator before surgery.

9. HANDLING


To avoid scratching or damaging the implants, these should be handled with the utmost care by qualified personnel and in an environment where conditions of hygiene are controlled. The implants should be kept in their undamaged packages.

10. SURGICAL TECHNIQUE

The surgeon should be fully familiar with the surgical technique. Supplementary information about the surgical techniques (brochure and video) and products are available on request. Careful preoperative planning is essential.

MECTALOCK PEEK & MECTALOCK C ANCHOR (Size 2.4 & 2.9mm)

- Pass the provided suture through soft tissue with surgeon's preferred suture-passing technique; let the two strands of the suture come out from the patient through the arthroscopic portal
- Find the desired location for the placement of the MectaLock PEEK or C anchor; insert a percutaneous guidewire (max Ø1.5mm) through the proper arthroscopic portal to reach the desired location
- A dedicated aimer (slotted or non-slotted) can be used as a drill guide and to aid in pointing the location of the pilot hole in the bone. Insert the Medacta Arthroscopic Aimer coupled with a cannulated obturator over the guidewire; the surgeon can choose between different aimer's tip according to his preference and the anatomical structures to be treated. Aimers and obturators are provided in two different lengths:


Aimers/obturators and drills can be coupled according to the following table:

AIMER		OBTURATOR		REUSABLE DRILL BIT	
REFERENCE	DESCRIPTION	REFERENCE	DESCRIPTION	REFERENCE	DESCRIPTION
05.11.10.0030	Crown Tip Long	05.11.10.0033	Obturator Long	05.11.10.0035 05.11.10.0036	Drill Ø2.4 Long Drill Ø2.9 Long
05.11.10.0052	Crown Tip Long Slotted				
05.11.10.0031	Fishmouth Tip Long				
05.11.10.0053	Fishmouth Tip Long Slotted				
05.11.10.0032	Fork Tip Long	05.10.10.0013	Obturator Short	05.10.10.0015 05.10.10.0016	Drill Ø2.4 Short Drill Ø2.9 Short
05.11.10.0054	Fork Tip Long Slotted				
05.10.10.0010	Crown Tip Short				
05.10.10.0024	Crown Tip Short Slotted				
05.10.10.0011	Fishmouth Tip Short				
05.10.10.0025	Fishmouth Tip Short Slotted				
05.10.10.0012	Fork Tip Short	05.10.10.0026	Fork Tip Short Slotted		
05.10.10.0026	Fork Tip Short Slotted				

- Using the corresponding compatible reusable drill bit, create a hole in the bone in the desired location for the placement of the MectaLock PEEK or C anchor
- Make sure to drill until the distal depth mark of the drill is aligned with the distal depth mark of the aimer or until the drill is mechanically blocked against the aimer


Load the MectaLock PEEK or C anchor with the sutures coming out of the portal: while firmly holding the two strands of the suture with one hand, align them with the slot of the anchor's tip. Gently push the anchor against the suture to let it clip into the tip's eyelet.

Rotate the anchor clockwise by 90°


Align the main axis of the driver with the suture (1) and then, while still holding the suture gently tensioned, slide the anchor over the suture (2) towards the patient through the arthroscopic portal, until the pilot hole is reached.

Once reaching the pilot hole, regulate the suture tension according to the surgeon's preference and tap the plastic handle of the driver until the MectaLock PEEK or C anchor is correctly placed. The correct depth is achieved when the distal laser marking is flush with the bone surface.


- Unscrew the disposable driver
- Palpate the repaired soft tissue using an arthroscopic probe to ensure that the repair is secure
- Cut and remove sutures in excess

MECTALOCK PEEK & MECTALOCK C ANCHOR (Size 4.5, 5.5 & 6.5 mm)

Using the corresponding compatible reusable punch, create a socket in the bone in the desired location for the placement of the MectaLock PEEK or MectaLock C Suture anchor.

- Make sure to punch until the distal depth mark of the punch is flush with cortical bone.


- Load the MectaLock PEEK or C anchor with the sutures coming out of the portal: while firmly holding the two strands of the suture with one hand, align them with the slot of the anchor's tip.
- Gently push the anchor against the suture to let it clip into the tip's eyelet.
- Rotate the anchor clockwise by 90°.

- Align the main axis of the driver with the suture (1) and then, while still holding the suture gently tensioned, slide the anchor over the suture (2) towards the patient through the arthroscopic portal, until the socket hole is reached.
- Once reaching the socket hole, regulate the suture tension according to the surgeon's preference and tap the plastic handle of the driver until the MectaLock PEEK or C anchor is correctly placed. The correct depth is achieved when the distal laser marking is flush with the bone surface.
- Unscrew the disposable driver.
- Palpate the repaired soft tissue using an arthroscopic probe to ensure that the repair is secure.
- Cut and remove sutures in excess.


MECTALOCK TI SUTURE ANCHOR

- Find the desired location for the placement of the MectaLock TI anchor; insert the preloaded driver with the MectaLock TI anchor through the proper arthroscopic portal to reach the desired location.
- Gently tap on the back handle of the disposable driver to let the MectaLock TI anchor's tip pierce the bone in the chosen spot; while keeping a gentle axial force on the disposable driver, screw the MectaLock TI anchor into the bone until the distal laser marking on the driver is flush with the bone surface.

Unlock the sutures from the plastic handle and gently disengage the driver's tip from the MectaLock TI anchor, pulling the driver out through the arthroscopic portal.

- Manage the 2x USP#2/EP#5 HS Fiber 39" UHMWPE sutures through soft tissue with surgeon's preferred suture-passing and knotting technique
- Palpate the repaired soft tissue using an arthroscopic probe to ensure that the repair is secure
- Cut and remove suture in excess

11. POSTOPERATIVE CARE AND FOLLOW-UP

The surgeon should caution the patients to control their level of activity and avoid excessive loads on the operated limb. Moreover, the surgeon should make the patients aware of the precautions to be taken in terms of exercises, treatments and limitations on activities, any limitations reported on the label, as well as exposure to magnetic fields. The patient must be told that implants can affect the results of computer tomography (CT) or magnetic resonance imaging (MRI) scans. Periodic follow-up are recommended to make comparisons with the immediate postoperative condition and anticipate implant related complications. Excessive physical activity, and operated limb traumas may cause early failure of the implant. If the case occurs, it is necessary to place the patient under supervision, evaluate the possible progression of the deterioration, and weigh the benefit of early revision.

12. ADVERSE EVENTS AND COMPLICATIONS

Adverse events that can occur in the reattachment of soft tissue to bone in orthopedic surgical procedures:

- Infection, both deep and superficial
- Allergies, mild inflammatory and foreign body reactions to implant material

Some adverse effects can ultimately lead to death.

General complications include:

- Venous thrombosis with/without pulmonary embolism
- Cardiovascular or pulmonary disturbances
- Haematomas
- Systemic allergic reactions
- Systemic pain

13. PACKAGING

The MectaLock Suture Anchors are supplied in single-use packages. For components delivered sterile, the sterilization method is indicated on the label. The expiration date must be checked on the label as well as the package integrity to ensure that sterility of the contents has not been compromised. If the package is damaged or has been previously opened, do not use the component. Do not re-sterilize.

14. STORAGE

The packages must be stored in a cool, dry place, away from light.

For MectaLock C Suture anchor (Composite):

The packages must be stored in a cool, dry place, away from light and at a maximum temperature of 25°C. If the product has been exposed to excessively high temperatures, the square on the temperature indication label on the outer packaging changes color from grey to black. MectaLock C Suture anchor with a black indicator must not be used.

15. TRADEMARKS

Medacta is registered trademark of Medacta International SA, Castel San Pietro, Switzerland.

The reference text is the English text.

This document is not intended for the US market.

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