## INTERFERENCE SCREW: MECTASCREW

### CAUTION

For a U.S. specific package insert, please go to [www.medacta.com](http://www.medacta.com) or call toll free 800-901-7836.

#### SYMBOLS / SIMBOLI / PICTOGRAMMES / PIKTOGRAMME / PICTOGRAMAS / PICTOGRAMAS / PICTOGRAMMES / ΕΙΚΟΝΟΣΥΜΒΟΛΑ

	Do not reuse / Non riutilizzare / Niet hergebruiken / Nicht wiederverwenden / No reutilizable / Não reutilizar / Ne pas réutiliser / Μην επαναχρησιμοποιείτε
	Do not re-sterilize / Non risterilizzare / Niet re-steriliseren / Bitte nicht re-sterilisieren / No re-sterilize / Não re-sterilizar / Ne pas restériliser / Να μην αποστειρωθεί ξανά
	Caution, read the accompanying documents / Attenzione, consultare i documenti di accompagnamento / Opgelet, raadpleeg de bijgeleverde documenten / Achtung, Begleitdokumente beachten / Atención, consulte los documentos que se acompañan / Atenção, consultar os documentos de acompanhamento / Attention, consulter les documents d'accompagnement / Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα
	Consult instructions for use / Consultare le istruzioni di utilizzo / Gebruiksaanwijzing / Bitte lesen Sie die Gebrauchsanweisung / Consulte las Instrucciones antes de su uso / Consulte as instruções de uso / Consultez les instructions d'utilisation / Οδηγίες χρήσης
	Do not expose to sunlight / Non esporre alla luce del sole / Niet blootstellen aan zonlicht / Vor Sonnenlicht schützen / No exponer a la luz solar / Não expor à luz solar / Ne pas exposer à la lumière du soleil / Μην εκθέτετε σε ηλιακό φως
	Store in a dry place / Conservare in luogo asciutto / Droog bewaren / Trocken aufbewahren / Conservar en lugar seco / Conservar em local seco / Conserver au sec / Φυλάσσετε σε ξηρό χώρο
	Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Niet gebruiken als verpakking beschadigd is / Bitte nicht verwenden, falls Verpackung beschädigt ist / No utilize si el envoltorio está dañado / Não utilizar se a embalagem estiver violada / Ne pas utiliser si l'emballage est abîmé / Να μην χρησιμοποιηθεί εάν η συσκευασία είναι κατεστραμμένη
	Use by / Utilizzare entro il / Te gebruiken voor / Zu verwenden bis / Utilizar antes de / A utilizar antes de / A utiliser avant / Ημερομηνία λήξης
	Lot number / Numero di lotto / Nummer van de partij / Waren-Lot-Nummer / Número de lote / Número de lote / Numéro de lot / Αριθμός παρτίδας
	Reference number / Riferimento commerciale / Handelsreferentie / Referenznummer / Numero de Referencia / Referência comercial / Référence commerciale / Νούμερο παραπομπής
	Sterilized by irradiation / Sterilizzato mediante irradiazione / Gesteriliseerd door doorstraling / Sterilisiert durch Bestrahlung / Esterilizado por irradiación / Esterilizado por irradiação / Stérilisé par irradiation / Αποστειρωμένο με ακτινοβολία
	Manufacturer / Produttore / Fabrikant / Hersteller / Fabricante / Fabricante / Fabricant / Κατασκευαστής
	RxOnly / Solo prescrizione medica / RX Enkel / RX Only / Sólo RX / Somente RX / Sur prescription seulement / RX μόνον
	MR Safe / Sicuro per RM / MR veilig / MR geeignet / Seguro para RM / Seguro para RM / Compatible avec IRM / Ασφαλές στην μαγνητική τομογραφία
	Temperature limit
	Double sterile barrier system / Sistema a Doppia Barriera Sterile / Dubbele steriele verpakking / Doppeltes steriles Barriersystem / Sistema de doble barrera estéril / Sistema de Dupla Barreira Estéril / Système à barrière stérile double / Σύστημα διπλού στείρου εμποδίου

## Index

ENGLISH: MECTASCREW INTERFERENCE SCREW - INSTRUCTIONS FOR USE .....	3
ITALIANO: VITE DI INTERFERENZA MECTASCREW - ISTRUZIONI PER L'USO.....	6
NEDERLANDS: MECTASCREW INTERFERENTIESCHROEF- GEBRUIKSINSTRUCTIES .....	9
DEUTSCH: MECTASCREW INTERFERENZSCHRAUBE - GEBRAUCHSANWEISUNGEN.....	12
ESPAÑOL: TORNILLO INTERFERENCIAL MECTASCREW - INSTRUCCIONES DE USO .....	15
PORTUGUÊS: PARAFUSO DE INTERFERÊNCIA MECTASCREW - INSTRUÇÕES DE UTILIZAÇÃO.....	18
FRANÇAIS: VIS D'INTERFÉRENCE MECTASCREW - INSTRUCTIONS D'UTILISATION .....	21
ΕΛΛΗΝΙΚΑ: ΒΙΔΑ ΣΤΕΡΕΩΣΗΣ MECTASCREW - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ.....	24

## ENGLISH: MECTASCREW INTERFERENCE SCREW - INSTRUCTIONS FOR USE

Important notice: the device(s) can be prescribed and implanted only by a doctor legally authorized to perform this type of surgery.

### 1. GENERAL

Before any surgery, the surgeon must be familiar with the product literature and operative technique and must read carefully these instructions for use. Patient selection is as important as implant placement or positioning. Unsuitable functional requirements may contribute to reduce the implant life. The warnings must be heeded, and the instructions for use must be strictly followed.

### 2. PRODUCT DESCRIPTION

These instructions for use are intended for all Medacta MectaScrew Interference Screw implants. Medacta MectaScrew Interference Screw has been developed to provide interference between a ligament graft and a bone tunnel in the ligament reconstruction surgery.

DEVICE NAME	MATERIAL CODE	MATERIAL TRANSCRIPTION
MectaScrew Interference Screw PEEK	PEEK	Polyether Ether Ketone
MectaScrew B Interference Screw	PLDLLA	Poly(L-lactide-co-D,L-lactide)
MectaScrew C Interference Screw	COMPOSITE	Poly(L-lactide-co-D,L-lactide) + $\beta$ -TCP

### 3. INTENDED USE / INDICATIONS

The purpose of the MectaScrew Interference Screw is the reconstructive treatment of ruptured anterior and posterior cruciate ligaments by means of an auto- or allograft.

### 4. CONTRAINDICATIONS

- Degenerative osteopathies
- Local Bone tumors
- Osteomata in the area in which the implant has to be placed
- Deformities of the bone, or general conditions of the bone, which preclude the implantation in the opinion of a physician
- Systemic diseases and metabolic disorders that may compromise the outcome of the surgery

### 5. WARNINGS AND PRECAUTIONS

The success of the operation depends on compliance with the operative technique supplied, and the proper use of the dedicated instruments specially designed for that range of implants. Malpositioning may reduce implant longevity and lead to early implant failure. One Medacta label for each device must be properly applied to the Medacta Implant Card provided and delivered to the patient.

Ensure Medacta MectaScrew Interference Screws are used with the special screwdriver designed for this purpose. The screwdriver must be inserted into the screw as far as it will go in order to prevent damage to the screw and its head during insertion.

#### MRI Compatibility

The Medacta MectaScrew Interference Screw is MR Safe.

## 6. MEDACTA INTERNATIONAL IMPLANTS

Medacta® International is not responsible for the use of its implant components in combination with a component from another manufacturer (unless otherwise specified by Medacta® International in the surgical technique), therefore we advise against such a use. The components of a MectaScrew Interference Screw should never be reimplanted. While an implant may appear undamaged, microscopic imperfections may occur and cause implant failure. The operating surgeon must be aware that even a very small superficial damage, caused for instance by a sharp tool or electrocautery, can have an influence on the endurance of the device and can lead to fracture. Medacta MectaScrew Interference Screw products are non-pyrogenic.

## 7. RISK FACTORS

The following risk factors, individually or together, may result in poor clinical outcomes:

- Inadequate bone quality, (e.g. osteoporosis, previous cruciate ligament operation)
- Systemic diseases or metabolic disorders
- History of infections or recurrent falls
- Osteoporosis or osteomalacia
- Drug dependence and abuse of alcohol and medicaments
- Mental incapacity of patient to understand the instructions of the physician and to comply with them
- Adiposity or patient overweight, leading to excessive strain on the fixation

## 8. INSTRUCTIONS FOR USE

### PREOPERATIVE PHASE

The surgeon should verify possible patient physical limitations and mental deficiencies and he should also discuss with the patient all the details of the procedure and implant. The discussion should consider the limitations of the procedure and the constraints imposed by the selected implant. The factors which could limit the performance and stability of the implant, e.g. level of activity, should be set out to improve the patient's chances to avoid complications. The necessity to follow the postoperative instructions given by the surgeon should be fully understood by the patient. A stock of sterile implants of suitable sizes must be available and checked by the operator before surgery.

## 9. HANDLING

To avoid scratching or damaging the implants, these should be handled with the utmost care by qualified personnel and in an environment where conditions of hygiene are controlled. The implants should be kept in their undamaged packages.

## 10. SURGICAL TECHNIQUE

The surgeon should be fully familiar with the surgical technique. Supplementary information about the surgical techniques (brochure and video) and products are available on request. Careful preoperative planning, carried out on X-rays (and/or computer tomography (CT) or magnetic resonance imaging (MRI) scans), is essential.

## 11. POSTOPERATIVE CARE AND FOLLOW-UP

The surgeon should caution the patients to control their level of activity and avoid excessive loads on the operated limb. Moreover, the surgeon should make the patients aware of the precautions to be taken in terms of exercises, treatments and limitations on activities, any limitations reported on the label.

## 12. ADVERSE EVENTS AND COMPLICATIONS

Adverse events that can occur in reconstructive treatment of ligament ruptures using Medacta MectaScrew Interference Screw include:

- Infection, both deep and superficial;
- Allergies, mild inflammatory and foreign body reactions to implant material.

Some adverse events can ultimately lead to death

General complications include:

- Venous thrombosis with/without pulmonary embolism;
- Cardiovascular or pulmonary disturbances;
- Haematomas;
- Systemic allergic reactions;
- Systemic pain.

In addition, in case of use of bio-absorbable material (valid for MectaScrew B and MectaScrew C Interference Screw), adverse events that can occur include:

- In rare cases light and temporary tissue reaction may occur after surgery resulting in a joint effusion which regresses in a short period of time
- In very rare cases, slight expansion of the bone may occur around the implant site (osteolysis), potentially resulting in implant migration and/or loss of graft fixation

## 13. PACKAGING

The components of all MectaScrew Interference Screw are supplied in single-use packages. For components delivered sterile, the sterilization method is indicated on the label. The expiration date must be checked on the label as well as the package integrity to ensure that sterility of the contents has not been compromised. If the package is damaged or has been previously opened, do not use the component. Do not resterilize.

## 14. STORAGE

### 14.1 MECTASCREW PEEK INTERFERENCE SCREW

The packages must be stored in a cool, dry place, away from light.

### 14.2 MECTASCREW B AND MECTASCREW C INTERFERENCE SCREW

The packages must be stored in a cool, dry place, away from light and at a maximum temperature of 25°C. If the product has been exposed to excessively high temperatures, the square on the temperature indication label on the outer packaging changes color from grey to black. MectaScrew B and MectaScrew C Interference Screw with a black indicator must not be used.

## 15. TRADEMARKS

Medacta is registered trademark of Medacta International SA, Castel San Pietro, Switzerland.

The reference text is the English text

This document is not intended for the US market.

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