Distributed by:
Medacta International SA
Strada Regina CH- 6874 Castel San Pietro - Switzerland
Phone: +41 91 6966060 Fax: +41 91 6966066
info@medacta.ch - www.medacta.com

V.Mueller™ Sterile Distraction Screws
Manufactured by:
CareFusion
1500 Waukegan Rd - McGaw Park, IL 60085
Phone : 800-323-9088 - www.carefusion.com

All trademarks and registered trademarks are the property of their respective owners.

Mecta-C Cervical Distraction Sterile Pins


SYMBOLS / SIMBOLI / PICTOGRAMMES / PIKTOGRAMME / PICTOGRAMAS / PICTOGRAMMES

 Sterilized by irradiation / Sterilizzato mediante irradiazione / Gesteriliseerd door doorstraling / Sterilisiert durch Bestrahlung / Esterilizado por irradiación / Esterilizado por irradiação / Stérilisé par irradiation

 USA: solo Rx / Verenigde Staten: Enkel RX / USA. Rx Only / USA.: sólo Rx / EUA: só Raios-X / USA:seulement Rx

 Do not re-sterilize / Non risterilizzare / Niet resteriliseren / Bitte nicht re-sterilisieren / No re-sterilize / Não risterilizar / Ne pas restériliser

 Manufacturer / Produttore / Fabrikant / Hersteller / Fabricante / Fabricante / Fabricant / Κατασκευαστής:
CareFusion - 1500 Waukegan Rd - McGaw Park, IL 60085

 Care Fusion France 309 S.A.S. - 8 bis rue de la Renaissance, 44110 Chateaubriant - France
36-10389 ● 2013-06

 Product and packaging are not manufactured with natural rubber latex / Laex: Il Prodotto e il confezionamento non sono prodotti in Lattice naturale / Latex: Product en verpakking worden vervaardigd zonder natuurlijk rubber / Produkt und Verpackung werden nicht mit Naturkautschuklatex hergestellt / Látex: El producto y el embalaje no son fabricados en látex natural / Latex: O Produto e o acondicionamento não são produzidos em latex natural / Le produit et le conditionnement ne sont pas produit en latex naturel

ENGLISH Mecta-C RETRACTOR STERILE PINS - INSTRUCTIONS FOR USE

CATALOG NUMBERS

S-0088, S-0089, S-0090, S-0091

INDICATIONS FOR USE

Distraction screws are used in conjunction with orthopedic retractors (distractors) for retraction of vertebral bodies during surgical procedures.

HOW SUPPLIED

Distraction screws are supplied as sterile. Do not use if single use sterile package is damaged.

LIMITATIONS ON REPROCESSING

Distraction Screws are for single use only. Reuse may result in a non functional product or contribute to cross contamination.

WARNINGS

Devices shall be used in accordance with these Instructions for Use. Read all sections of this IFU prior to use. Improper use of these devices may cause serious injury. In addition, improper storage of the devices may render the devices non-sterile prior to patient use and may cause serious injury to the health care provider or the patient.

CAUTIONS

Use of a device for a task other than what it is intended for will usually result in a damaged or broken device. Prior to use, inspect devices to ensure proper function and condition. Do not use devices if they do not satisfactorily perform their intended function or if they have physical damage.

INSPECTION

Distraction screws with dull or damaged tips (dull, broken, bent, burrs, rough surfaces), or tarnished surfaces should not be used, but should be replaced immediately.

STORAGE

Devices should remain in sterilization packaging and be stored in a clean, dry environment. Dispose of distraction screws in accordance with accepted medical practice and all applicable local, state and federal laws and regulations.

WARRANTY

CareFusion offers a lifetime guarantee on every surgical device bearing the V. Mueller brand name (unless otherwise noted). Surgical devices are guaranteed to be free of functional defects in workmanship and materials when used normally for its intended surgical purpose. Any V. Mueller device proving to be defective will be replaced at no cost to the customer.

The reference text is the English text (from the manufacturer).

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