Publication Proof

Company Riverpoint Medical
 Customer _____

Number _____

Description Brochure FU

Title _____

Riverpoint Medical Approval

Each department must initial and date.

Art Department _____

Sales Department _____

Q.A. R.A. _____

Document Coordinator and DRF stamps

Customer Approval

The document is attached to this form. Examine carefully. Mark all changes and/or typographical errors. If you approve errors on this proof, Riverpoint Medical cannot be held responsible. Examine the entire job once again to ensure everything is correct.

If approving with changes, mark this box

Signature _____ Date _____

Customer signature represents approval of publication as shown. Required for private label only.

EN

PowerSuture™

Non-Absorbable surgical suture

PowerSuture Braided Sutures Instructions For Use

DESCRIPTION: PowerSuture sutures are non-absorbable, sterile, surgical sutures composed of ultra high molecular weight polyethylene (UHMW PE). These braided sutures are available uncoated or coated uniformly with polybutylene adipate or silicone to increase surface lubricity, thereby enhancing the handling characteristics, ease of passage through tissue, and knot run-down properties for security.

PowerSuture sutures are inert and elicit only minimal local tissue reaction. The suture braid is available undyed (white), dyed blue, dyed black or with trace filaments of black, blue or green suture for color.

PowerSuture sutures meet all requirements established by the United States Pharmacopeia (USP) for non-absorbable surgical sutures except where indicated on the individual unit package.

INDICATIONS: PowerSuture sutures are indicated for use in general soft tissue approximation and/or ligation and the use of allograft tissues for orthopedic surgeries, but are not for use in cardiovascular procedures. PowerSuture sutures are intended for one-time use only, and are not to be re-sterilized.

ACTIONS: PowerSuture sutures elicit a minimal acute inflammatory reaction in tissues, followed by gradual encapsulation of the suture by fibrous connective tissue. PowerSuture sutures are not absorbed, nor is there any significant change in tensile strength retention known to occur in vivo.

CONTRAINDICATIONS: None known.

WARNINGS: As with any foreign body, prolonged contact of any suture with salt solutions may result in calculus formation. Discard opened or unused sutures. Do not resterilize. Users should be familiar with surgical procedures and techniques involving sutures before employing PowerSuture non-absorbable sutures for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

PRECAUTIONS: In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments.

ADVERSE REACTIONS: Adverse effects associated with the use of this device include: wound dehiscence, calculi formation when prolonged contact with salt solutions occur, enhanced bacterial infectivity, minimal acute inflammatory reaction, and transitory local irritation.

HOW SUPPLIED: PowerSuture sutures are available undyed (white), dyed blue, dyed black, or with trace filaments composed of polyester, nylon, or polypropylene, dyed black, blue, or green to add visual color to the suture. PowerSuture is available coated or uncoated, in sizes 7 (9 metric) through 6-0 (0.7 Metric). PowerSuture Tape is a flat braid PowerSuture suture offered in widths from 1MM to 5mm. The sutures are supplied sterile, in pre-cut lengths, non-needed or needed. PowerSuture sutures are available in up to one dozen box quantities. This device was designed, tested and manufactured for single patient use only. Reuse or reprocessing of this device may lead to its failure and subsequent patient injury. Reprocessing and/or resterilization of this device may create the risk of contamination and patient infection. Do not reuse, reprocess or resterilize this device.

CAUTION: Federal (USA) law restricts this device to sale and use by, or on the order of, a physician or veterinarian.

CAUTION: Federal (USA) law restricts this device to sale and use by, or on the order of, a physician or veterinarian.

CAUTION: Federal (USA) law restricts this device to sale and use by, or on the order of, a physician or veterinarian.

DO NOT RESTERILIZE ATTENTION READ INSTRUCTIONS

ETHYLENE OXIDE STERILIZATION SINGLE USE ONLY

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