M.U.S.T. SACRAL ILIAC SCREWS SYSTEM

CAUTION

For a U.S. specific package insert, please go to www.medacta.com or call toll free 800-901-7836.

SYMBOLS / SIMBOLI / PICTOGRAMMES / PIKTOGRAMME / PICTOGRAMAS / PICTOGRAMAS / PICTOGRAMMES

Do not reuse / Non riutilizzare / Niet hergebruiken / Nicht wiederverwenden / No reutilizable / Não reutilizar / Ne pas réutiliser

Do not re-sterilize / Non risterilizzare / Niet re-steriliseren / Bitte nicht re-sterilisieren / No re-sterilize / Não reesterilizar / Ne pas restériliser

Caution, read the accompanying documents / Attenzione, consultare i documenti di accompagnamento / Opgelet, raadpleeg de bijgeleverde documenten / Achtung, Begleitdokumente beachten / Atención, consulte los documentos que se acompañan / Atenção, consultar os documentos de acompanhamento / Attention, consulter les documents d'accompagnement

Consult instructions for use / Consultare le istruzioni di utilizzo / Gebruiksaanwijzing / Bitte lesen Sie die Gebrauchsanweisung / Consulte las Instrucciones antes de su uso / Consulte as instruções de uso / Consultez les instructions d'utilisation

Do not expose to sunlight / Non esporre alla luce del sole / Niet blootstellen aan zonlicht / Vor Sonnenlicht schützen / No exponer a la luz solar / Não expor à luz solar / Ne pas exposer à la lumière du soleil

Store in a dry place / Conservare in luogo asciutto / Droog bewaren / Trocken aufbewahren / Conservar en lugar seco / Conservar em local seco / Conserver au sec

Do not use if package is damaged / Non utilizzare se la confezione è danneggiata / Niet gebruiken als verpakking beschadigd is / Bitte nicht verwenden, falls Verpackung beschädigt ist / No utilize si el envoltorio está dañado / Não utilizar se a embalagem estiver violada / Ne pas utiliser si l'emballage est abimé

Use by / Utilizzare entro il / Te gebruiken voor / Zu verwenden bis / Utilizar antes de / A utilizar antes de / A utiliser avant

Lot number / Numero di lotto / Nummer van de partij / Waren-Lot-Nummer / Número de lote / Número de lote

Reference number / Riferimento commerciale / Handelsreferentie / Referenznummer / Numero de Referencia / Referência comercial / Référence commerciale

Sterilized with ethylene oxide / Sterilizzato con ossido di etilene / Gesteriliseerd met ethyleenoxide / Ethylenoxid-sterilisiert / Esterilizado con óxido de etileno / Esterilizado com óxido de etileno / Stérilisé à l'oxyde d'éthylène

Sterilized by irradiation / Sterilizzato mediante irradiazione / Gesteriliseerd door doorstraling / Sterilisiert durch Bestrahlung / Esterilizado por irradiación / Esterilizado por irradiação / Stérilisé par irradiation

Manufacturer / Produttore / Fabrikant / Hersteller / Fabricante / Fabricante / Fabricant

MR Conditional

Index

ENGLISH: M.U.S.T. SACRAL ILIAC SCREWS SYSTEM - INSTRUCTIONS FOR USE.....	3
ITALIANO: SISTEMA DI VITI PER FISSAZIONE SACRO-ILIACA M.U.S.T. - ISTRUZIONI PER L'USO.....	8
NEDERLANDS: M.U.S.T. SCHROEFSTEEM BIJ ILIO-SACRALE AANDOENINGEN - GEBRUIKSINSTRUCTIES	14
DEUTSCH: M.U.S.T. SAKROILIAKALGELENKSYSTEM - GEBRAUCHSANWEISUNG	20
ESPAÑOL: SISTEMA DE TORNILLO SACROILÍACO M.U.S.T. - INSTRUCCIONES DE USO	26
PORTUGUÊS: M.U.S.T. SISTEMA DE PARAFUSOS SACROILÍACOS - INSTRUÇÕES DE UTILIZAÇÃO	32
FRANÇAIS : SYSTÈME DE VIS M.U.S.T. SACRO-ILIAQUES - MODE D'EMPLOI.....	38

ENGLISH: M.U.S.T. SACRAL ILIAC SCREWS SYSTEM - INSTRUCTIONS FOR USE

Important notice: the device(s) can be prescribed and implanted only by a doctor legally authorized to perform this type of surgery.

1. GENERAL

Before any surgery, the surgeon M.U.S.T. be familiar with the sales literature and operative technique and M.U.S.T. carefully read these instructions for use. Patient selection is as important as implant placement and positioning. Obesity or unsuitable functional requirements may generate exceptional stress and reduce the implant life. The warnings M.U.S.T. be heeded and the instructions for use M.U.S.T. be strictly followed.

2. INTENDED USE

The M.U.S.T. Sacro Iliac Implant System [M.U.S.T. SI] is designed for the sacroiliac joint fusion.

3. PRODUCT DESCRIPTION

The M.U.S.T. Sacro Iliac screw system includes:

- Titanium screws with head (Standard screws) or without (Headless screws) plasma spray coated with porous and rough Hydroxyapatite.
- Titanium washers to be used with head screw (Standard screw)

4. INDICATIONS FOR USE

The M.U.S.T. SI Implant System is intended for sacroiliac joint fusion for skeletally mature patients suffering from sacroiliac joint disruptions, degenerative sacroiliitis and degenerative sacro-iliac arthritis, secondary to pelvic disruption.

5. CONTRAINDICATIONS

- Deformities or anatomic variations that prevent or interfere with implant placement.
- Bone tumor involving the site of operation.
- Active infection at treatment site.
- Intolerance / Allergy to the materials used in the manufacture of this device.
- Any active or suspected latent infection or marked local inflammation in or about the affected area.
- Compromised vascularity that would inhibit adequate blood supply to the operative site.
- Patients with, fever, tumors, elevated white blood count, mental illness and other medical conditions which would prohibit beneficial surgical outcome.
- Patients having inadequate tissue coverage over the operative site or inadequate bone stock or quality that cannot provide adequate support and/or fixation of the devices.
- Implant utilization that would interfere with anatomical structures or physiological performance.
- Any neuromuscular disorder which could create an unacceptable risk of fixation failure or complications in post-operative care.
- Other medical or surgical conditions which would preclude the potential benefit of surgery.
- Reuse or multiple uses.
- Unstable fracture of sacrum and or ilium involving the sacroiliac joint.

6. WARNINGS

- A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where other patient conditions may compromise the results.
- Preoperative and operating procedures, including knowledge of the surgical techniques and correct selection and placement of the implants are important considerations in the successful utilization of the system by the surgeon. Furthermore, the proper selection and the compliance of the patient will greatly affect the results. Patients who smoke have been shown to have a reduced incidence of bone fusion. These patients should be advised of this

fact and warned of this consequence. Obese, malnourished and /or alcohol /drug abuse patients and those with poor muscle and bone quality and /or nerve paralysis are also poor candidates for spinal fixation.

- The following circumstances may reduce the chance of a successful outcome: signs of local inflammation; fever or leucocytosis; metal/polymer sensitivity/allergies to the implant material; grossly distorted anatomy due to congenital abnormalities; osteopenia and/or osteoporosis; unsuitable or insufficient bone support.
- CORRECT SELECTION OF THE IMPLANT IS EXTREMELY IMPORTANT. The potential for satisfactory fixation is increased by the selection of the proper size, shape, and design of the implant. While proper selection can help minimize risks, the size and shape of human bones present limitations on the size, shape and strength of implants.
- The M.U.S.T. SI screw system implants should never be re-implanted. While an implant may appear undamaged, microscopic imperfections may occur and cause implant failure.

One Medacta label for each device M.U.S.T. be properly applied to the Medacta Implant Card provided, and delivered to the patient.

MRI Compatibility

Non-clinical testing has demonstrated the M.U.S.T. SI Sacro Iliac Implant System is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3 T.
- Maximum spatial field gradient of 6,400 gauss/cm.
- Maximum Force Product 180,000,000 G²/cm (108 T²/m)
- The theoretically estimated maximum whole body average specific absorption rate (SAR) of <2 W/kg (Normal Operating Mode).

Under the scan conditions defined above, the M.U.S.T. SI Sacro Iliac Implant System is expected to produce a maximum temperature rise of less than:

- 2.4°C (2W/kg, 1.5Tesla) RF-related temperature increase with a background temperature increase of 0.8°C (2W/kg, 1.5Tesla)
- 2.1°C (2 W/kg, 3 Tesla) RF-related temperature increase with a background temperature increase of 0.2°C (2 W/kg, 3 Tesla) after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 69.3 mm from the M.U.S.T. SI Sacro Iliac Implant System when imaged with a gradient echo pulse sequence and a 3 T MRI system.

CAUTION

Patient safety is ensured up to a static magnetic field of 3T. No further analysis has been performed for higher values.

7. PRECAUTIONS

- The use of the M.U.S.T. SI screw system should only be undertaken after the surgeon has become thoroughly knowledgeable about anatomy and biomechanics, has had experience with procedures and has had hands-on training in the use of this device.
- The M.U.S.T. SI screw system should not be implanted in patients with severe osteoporosis or osteopenia.
- Correct handling of the implants is extremely important.
- Alterations will produce defects in surface finish and internal stresses which may become the focal point for eventual breakage of the implant. Bending of screws will significantly decrease the fatigue life and may cause failure.
- All the M.U.S.T. SI screw system devices supplied sterile should be handled with appropriate precautions to maintain sterility.
- The important medical information given in this document should be conveyed to the patient.
- An explanted metal implant should never be re-implanted. Even though the device appears undamaged, it may have small defects and internal stress patterns which may lead to early breakage.
- Pregnancy is not a contraindication for the implants here contemplated, but diagnostic and operative procedures on pregnant women should be undertaken after a comprehensive medical evaluation.
- Women of childbearing potential should be cautioned that vaginal delivery of a fetus may not be advisable following SI joint fusion. If pregnancy occurs, the woman should review delivery options with her obstetrician.
- Although not absolutely contraindicated, conditions to be considered as potential factors for not using this device include severe bone resorption, osteomalacia, and severe osteoporosis.

8. ADVERSE EFFECTS

Potential complications and adverse effects for this system are similar to those of other spinal instrumentation systems, and include, but are not limited to:

- Early or late loosening of any or all of the components.
- Foreign body allergic reaction or metal sensitivity to implants.
- Infection.
- Loss of neurological function including paralysis (partial or complete), radiculopathy, and/or the development or continuation of pain, numbness, spasms, or sensory loss.
- Cauda equina syndrome, neurological deficits, paraplegia, reflex deficits, irritation, and/or muscle loss.
- Scar formation possibly causing neurological compromise or compression around nerves and/or pain.
- Neurological injury/deficit which may range from paresthesias to muscle paralysis, loss of rectal or bladder sphincter control, radiculopathies.
- Fracture, micro-fracture, resorption, damage, or penetration of any spinal bone.
- Non-union (pseudarthrosis), delayed union, mal-union.
- Cessation of any potential growth of the operated portion of the spine.
- Cardiovascular system compromise.
- Decrease in bone density due to stress shielding.
- Device bending, disassembly, fracture, loosening, migration.
- Gastrointestinal complications (i.e. ileus or bowel perforation).
- Hemorrhage.
- Incisional complications (i.e. dehiscence, hematoma).
- Malfunction /Malposition of the fixation device.
- Organ, connective tissue or nerve damage.
- Osteoarthritis.
- Pain, discomfort or abnormal sensation due to device presence.
- Persistent low back pain.
- Reproductive system compromise.
- Screw back-out or breakage possibly leading to local pain, perforation or irritation of adjacent structures.
- Loss of or increase in joint mobility or function.
- Inability to perform the activities of daily living.
- Death.
- Potential difficulty in delivering fetus vaginally due to device-related restriction of SI joint stretching.

Adverse events may necessitate re-operation or revision. A revision is a procedure which adjusts or in any way modifies the original implant configuration, e.g. adjusting the position of the original configuration or removal of implants with their subsequent replacement. A removal is a procedure which removes one or more implants of the original device configuration without any replacement. A re-operation is a procedure which involves any surgical procedure at the involved sites which does not remove, modify or add any implants.

INSTRUCTIONS FOR USE

9. IMPLANT SELECTION:

- The selection of the proper size, shape and design of the implant for each patient is crucial to the success of the procedure. Surgical implants are subject to repeated stress during use and their strength is limited by the need to adapt the design to the human anatomy. Unless great care is taken in patient selection, placement of the implant and postoperative management to minimize stresses on the implant, such stresses may cause material fatigue and consequent breakage of the implant before the fusion process is complete, which may result in further injury or the need to remove the implant prematurely.
- Excessive loads, such as excessive torque, tensile or compression load applied to long handle insertion tools attached to the implant or direct application of loads to a small area of the devices can damage the implant/instrument interface.
- Once the fusion has healed, the surgeon and patient should carefully weigh the risks and benefits if considering the removal of the screw fixation system.
- Improper selection, placement, positioning and fixation of these devices may result in unusual stress conditions reducing the service life of the implant.
- The surgeon is to be thoroughly familiar with the surgical procedure, instruments and implant characteristics prior to performing surgery.

Periodic follow-up is recommended to monitor the position and state of the implants, as well as the condition of the adjoining bone.

10. IMPLANT FIXATION

Installation and positional adjustment of implants M.U.S.T. only be done with special ancillary instruments and equipment supplied and designated by Medacta. In the interests of patient safety it is therefore recommended that Medacta implants are not used with implants from any other source.

11. PRE-OPERATIVE

The surgeon and operating theatre personnel should wear sterile surgical gloves. Under no circumstances should the components come into contact with hard objects (unless otherwise specified in the surgical technique). Before use, each component M.U.S.T. be visually inspected for imperfections. Special surgical instruments are required for the surgery. To avoid damaging the implants and the instruments, these should be handled with the utmost care by qualified personnel and in an environment where conditions of hygiene are controlled. The correct usage of instruments is described in the surgical technique. Instruments M.U.S.T. be inspected visually before the operation. Distorted or damaged instruments may result in malposition of the implant or implant failure.

- Only patients that meet the criteria described in the indications should be selected.
- Patient conditions and / or predispositions such as those addressed in the aforementioned contraindications should be avoided.
- Care should be taken in the handling and storage of the implant(s). They should not be scratched or damaged.
- The surgeon should be familiar with the various implants before use and should personally verify that all implants are present before the surgery begins.
- The size of the implant(s) for the case should be determined prior to beginning the surgery. An adequate inventory of implant sizes should be available at the time of surgery, including sizes larger and smaller than those expected to be used.
- All components and instruments should be cleaned and sterilized before use. Additional sterile instruments should be available in case of an unexpected need.

12. INTRA-OPERATIVE

- Any instruction manuals should be carefully followed.
- At all times, extreme caution should be used around the nerve roots. Damage to the nerves will cause loss of neurological functions.
- Breakage, slippage or misuse of instruments or implants may cause injury to the patient or operative personnel.
- The implant surface should not be scratched or notched since such actions may reduce the functional strength of the implant.
- Utilize an imaging system to facilitate surgery.
- A guide wire should first be used, followed by a sharp tap. Caution: Be careful that the guide-wire, if used, is not inserted too deep, becomes bent, and/or breaks. Ensure that the guide-wire does not advance during tapping or screw insertion. Remove the guide-wire and make sure it is intact. Failure to do so may cause the guide wire or part of it to advance through the bone and into a location that may cause damage to underlying structures.
- Do not overtap or use a screw that is either too long or too large. Overtapping, using an incorrectly sized screw, or accidentally advancing the guidewire during tap or screw insertion, may cause nerve damage, hemorrhage, or the other possible adverse events listed elsewhere in this package insert.

13. POSTOPERATIVE

The physician's postoperative directions and warnings to the patient and the corresponding patient compliance are extremely important.

- Detailed instructions on the use and limitations of the implant should be given to the patient. The patient M.U.S.T. be warned that loosening and / or breakage of the implant(s) are complications which may occur as a result of early or excessive weight-bearing, mechanical vibrations, muscular activity or sudden jolts or shocks.

- The patient M.U.S.T. be warned that bending, loosening and/or breakage of the device(s) are complications which may occur as a result of excessive or early weight-bearing or muscular activity. The risk of bending, loosening, or breakage of a device during postoperative rehabilitation may be increased if the patient is active, or if the patient is debilitated or demented. The patient should be warned to avoid falls or sudden jolts.
- The patient should be advised not to smoke tobacco or utilize nicotine products, or to consume alcohol or non-steroidal or anti-inflammatory medications such as aspirin during the bone graft healing process.
- It is important that the immobilization of the union is established and confirmed by radiographic examination. If a non-union develops or if the implants loosen, migrate and / or break, the implant should be revised and / or removed immediately before serious injury occurs. By the mechanism of fatigue, these stresses can cause the eventual bending, loosening, or breakage of the device(s). It is important that immobilization of the surgical site be maintained until firm bony union is established and confirmed by roentgen examination. If a state of non-union persists or if the components loosen, bend, and/or break, the device(s) should be revised and/or removed immediately before serious injury occurs. The patient M.U.S.T. be adequately warned of these hazards and closely supervised to insure cooperation until bony union is confirmed.
- Any retrieved implants should be treated in such a manner that re-use in another surgical procedure is not possible.
- Adequately instruct the patient of the appropriated postoperative care. The patient's ability and willingness to follow instructions is one of the most important aspects of successful healing.
- Instructions to the patient to reduce stress on the implants are an equally important part of the attempt to avoid the occurrence of clinical problems that may accompany fixation failure.

14. PACKAGING

The packaging of every components M.U.S.T. be intact upon receipt. If a loaner or consignment system is used, all sets should be carefully checked for completeness and all components including instruments should be carefully checked to ensure that there is no damage prior to use. Damaged packages or products M.U.S.T. not be used, and should be returned to Medacta International or to the local distributor. M.U.S.T. SI screws are available in sterile packaging. For the implants provided sterile, the sterilization method is indicated on the label. The expiration date M.U.S.T. be checked on the label as well as the package integrity to ensure that sterility of the contents has not been compromised. If a sterile package is damaged or has been previously opened, do not use the implant.

15. INSTRUMENTS

Instruments are supplied non-sterile and M.U.S.T. be cleaned and sterilized prior to use. Recommended cleaning, decontamination and sterilization instructions are provided on www.medacta.com.

16. STORAGE

The packages M.U.S.T. be stored in a cool, dry place, not exposed to sunlight.

17. PRODUCT COMPLAINTS

Any healthcare professional (e.g. customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction in the product quality, identity, durability, reliability, safety, effectiveness and / or performance should notify Medacta or the distributor. Furthermore, if any of the implanted devices ever malfunctions (i.e. does not meet any of its performance specifications or otherwise does not perform as intended) or is suspected of doing so, the distributor should be notified immediately. If any Medacta product ever malfunctions and may have caused or contributed to the death or serious injury of a patient, the distributor should be notified immediately by telephone, fax or written correspondence. When filing a complaint, please provide the implant(s) name and number, lot number(s), your name and address, the nature of the complaint and notification of whether a written report from the distributor is requested.

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