


.U.S.T. SI

SACRO-ILIAC JOINT SCREW SYSTEM

A COMPLETE SYSTEM WITH DIFFERENT OPTIONS



Brochure


Joint

Spine

Sports Med

The M.U.S.T. SI, Medacta Universal Screw Technology Sacro Iliac System, is intended for sacroiliac joint fusion of skeletally mature patients suffering from sacroiliac **joint disruptions**, degenerative **sacroiliitis** and degenerative sacroiliac **arthritis**, secondary to pelvic disruption.

IMPLANT FEATURES



PRODUCT RANGE

Comprehensive range of implants to accommodate several patients anatomies.



M.U.S.T. SI Screws:

- 25-80mm length
- Ø8, 9, 10 mm

M.U.S.T. SI Washers:

- **Standard** Ø13, 15, 17, 20 mm. ROM 22°
- **Favoured angle** Ø15, 17, 20 mm. ROM up to 28°



swiss
made

