

FairFix

ADJUSTABLE BUTTON



ONE SIZE BUTTON WITH SIMPLIFIED SUTURE MANAGEMENT

EASY GRAFT PREPARATION AND SUTURE MANAGEMENT

ONE SIZE FITS ALL

3-POINT LOCKING DESIGN


**INTERCHANGEABLE FLIP
AND PULL SUTURE**

**COMPATIBLE WITH BUTTON
EXTENDER FOR REVISION CASES**

**GRAFT CRADLE TO PREVENT
SAWING EFFECT**

**STRONG CONSTRUCT WITH
MINIMAL DISPLACEMENT***

*Data on file



ORDERING INFORMATION

Order Number	Description
05.05.0091	FairFix Adjustable Button
05.05.0092	FairFix Adjustable Button Extender

RELATED PRODUCTS

Order Number	Description
05.05.10.0132	M-ARS ACL wires kit
05.05.10.0133	Ligament reconstruction wires kit
05.05.10.0135	Ø4.5mm cortical reamer pin
05.05.10.0136	Femoral Aimer AM Ø7-8 for cortical reamer pin
05.05.10.0137	Femoral Aimer AM Ø9-10 for cortical reamer pin
05.05.10.0145	Coated Ø2.4mm K-wire
05.05.10.0146	Coated Ø2.4mm passing wire

