

PowerKnot

HIGH STRENGTH SUTURE



HIGH STRENGTH SUTURE WITH INCREASED KNOT GRIP TECHNOLOGY


BETTER KNOT GRIP TECHNOLOGY

Special raw material treatment resulting in an optimized surface technology

MATERIAL

USP #2 made 100% out of Dyneema Purity (UHMWPE)

RUNNING DIRECTION INDICATION (RDI)



Available for black/blue/white and black/yellow/white color options

4 COLOR OPTIONS

- white/blue
- white/black
- black/blue/white
- black/yellow/white




All color options are available with or without needle

TECHNICAL DETAILS


- Suture length: 38" / 96.5 cm
- Needle specification: Taper Point (MO-6) / 26 mm / ½ circle
- PowerKnot High Strength Suture is provided single use sterile in a box of 12 each

"BETTER KNOT GRIP" PERFORMANCE

STATIC TENSILE STRENGTH



CYCLIC KNOT STABILITY



Data on file

- 4 knots (load @ 3mm)
- 4 knots (load @ final failure)
- 6 knots (load @ final failure)

ORDERING INFORMATION

Order Number	Description
05.03.0001	PowerKnot, USP #2, white/blue w/ needle
05.03.0002	PowerKnot, USP #2, white/blue
05.03.0003	PowerKnot, USP #2, white/black w/ needle
05.03.0004	PowerKnot, USP #2, white/black
05.03.0005	PowerKnot RDI, USP #2, black/blue/white w/ needle
05.03.0006	PowerKnot RDI, USP #2, black/blue/white
05.03.0007	PowerKnot RDI, USP #2, black/yellow/white w/ needle
05.03.0008	PowerKnot RDI, USP #2, black/yellow/white



SPORTSMED.MEDACTA.COM



Distributed by Medacta International SA | Strada Regina - 6874 Castel San Pietro - Switzerland | Phone +41 91 696 60 60 | Fax +41 91 696 60 66 | info@medacta.ch

Find your local dealer at: medacta.com/locations. All trademarks and registered trademarks are the property of their respective owners. This document is not intended for the US market. Please verify approval of the devices described in this document with your local Medacta representative.

ref: 99.123SM.1SALES-S rev. 00 Last update: September 2021