

MectaScrew PEEK
PEEK INTERFERENCE SCREW

MectaScrew C
COMPOSITE INTERFERENCE SCREW



STRESS-FREE INSERTION & ENHANCED GRAFT COMPRESSION

MATERIAL

2 different materials available:

- **MectaScrew PEEK**
- **MectaScrew C:**
Composite, Osteoconductive

70%	30%
PLDLLA	+ β -TCP

OPTIMIZED STRESS DISTRIBUTION

Torx connection for better load distribution

SIMPLIFIED AND SAFER INSERTION

- Blunt thread reducing graft laceration
- Sharp threaded tip design to facilitate screw insertion

DESIGN

- Special screwdriver tip design to simplify screw insertion
- Short, tapered screw design without compromising interference compression

SELF-RETAINING SCREWDRIVER TIP

A self-retaining screwdriver tip for the temporary securing of the screw, preventing accidental dropping

ORDERING INFORMATION

MECTASCREW PEEK INTERFERENCE SCREW			
Order Number	Description	Order Number	Description
05.05.0003	Interference Screw 6x15 T20	05.05.0013	Interference Screw 8x30 T25
05.05.0004	Interference Screw 6x20 T20	05.05.0014	Interference Screw 9x20 T25
05.05.0005	Interference Screw 6x25 T20	05.05.0015	Interference Screw 9x25 T25
05.05.0006	Interference Screw 7x15 T25	05.05.0016	Interference Screw 9x30 T25
05.05.0007	Interference Screw 7x20 T25	05.05.0017	Interference Screw 10x20 T25
05.05.0008	Interference Screw 7x25 T25	05.05.0018	Interference Screw 10x25 T25
05.05.0009	Interference Screw 7x30 T25	05.05.0019	Interference Screw 10x30 T25
05.05.0010	Interference Screw 8x15 T25	05.05.0020	Interference Screw 11x30 T40
05.05.0011	Interference Screw 8x20 T25	05.05.0021	Interference Screw 12x35 T40
05.05.0012	Interference Screw 8x25 T25		

MECTASCREW C INTERFERENCE SCREW			
Order Number	Description	Order Number	Description
05.05.0096	Composite Interference Screw 6x15 T20	05.05.0106	Composite Interference Screw 8x30 T25
05.05.0097	Composite Interference Screw 6x20 T20	05.05.0107	Composite Interference Screw 9x20 T25
05.05.0098	Composite Interference Screw 6x25 T20	05.05.0108	Composite Interference Screw 9x25 T25
05.05.0099	Composite Interference Screw 7x15 T25	05.05.0109	Composite Interference Screw 9x30 T25
05.05.0100	Composite Interference Screw 7x20 T25	05.05.0110	Composite Interference Screw 10x20 T25
05.05.0101	Composite Interference Screw 7x25 T25	05.05.0111	Composite Interference Screw 10x25 T25
05.05.0102	Composite Interference Screw 7x30 T25	05.05.0112	Composite Interference Screw 10x30 T25
05.05.0103	Composite Interference Screw 8x15 T25	05.05.0113	Composite Interference Screw 11x30 T40
05.05.0104	Composite Interference Screw 8x20 T25	05.05.0114	Composite Interference Screw 12x35 T40
05.05.0105	Composite Interference Screw 8x25 T25		

DEDICATED INSTRUMENTS			
Order Number	Description	Order Number	Description
05.05.10.0118	Cannulated Screwdriver Shaft T20	05.05.10.0122	Cannulated Screwdriver Shaft T40
05.05.10.0120	Cannulated Screwdriver Shaft T25	05.05.10.0124	Quick Connect Ratchet Handle cannulated

DISPOSABLES

Order Number	Description	Order Number	Description
05.05.10.0133	Ligament Reconstruction Wires Kit	05.05.10.0146	Coated Ø2.4mm Passing Wire
05.05.10.0145	Coated Ø2.4mm k-wire	05.05.10.0135	Ø4.5mm Cortical Reamer Pin

