

HIP ACCESS KIT



Product Catalogue

Joint

Spine

Sports Med

INDEX

1. HIP ACCESS KIT	3
1.1 Access Kit (Ref. 05.11.10.0041)	3
1.2 Reusable Handle for blades (Ref. 05.11.10.0040)	3

1. HIP ACCESS KIT

Medacta's Hip Access kit is designed to provide the surgeons with all the main instruments to gain the access to the surgical site. It is composed of two standard access needles, two nitinol guidewires and two arthroscopic blades. It is provide sterile and disposable. A reusable plastic handle for the blades is also available.



1.1 ACCESS KIT (Ref. 05.11.10.0041)

Access Needle

(2x)

Designed with a cannulation suitable to receive a $\varnothing 1.5\text{mm}$ guidewire and a working length of 180mm.

Provided with luer-lock connection and obturator.



Blades

Designed with two cutting edges to enlarge the access through the capsule's fibers and a special back-end to be connected with the reusable ergonomic handle.

Provided in two configurations: straight and curved.



Guidewire

(2x)

$\varnothing 1.5\text{mm}$, WL 430mm



1.2 REUSABLE HANDLE FOR BLADES (Ref. 05.11.10.0040)

Ergonomic handle with locking button to hold the Hip Access Kit blades.



Part numbers subject to change.

NOTE FOR STERILISATION

The disposable Access Kit is provide sterile. All other reusable instruments are not sterile upon delivery. They must be cleaned before use and sterilised in an autoclave noting the regulations of the country, EU directives where applicable and following the instructions for use of the autoclave manufacturer. For detailed instructions please refer to the document "Recommendations for cleaning decontamination and sterilization of Medacta International orthopaedic devices" available at www.medacta.com.



**REDEFINING BETTER
IN ORTHOPAEDICS
AND SPINE SURGERY**

MEDACTA.COM



Medacta International SA
Strada Regina - 6874 Castel San Pietro - Switzerland
Phone +41 91 696 60 60 - Fax +41 91 696 60 66
info@medacta.ch

Find your local dealer at: medacta.com/locations

All trademarks and registered trademarks are the property of their respective owners.
This document is intended for the US market.

Hip Access Kit
Product Catalogue

ref: 99.104SMH.180US
rev. 00

Last update: April 2019