

▲ MetaLoop

CONTINUOUS LOOP BUTTON



Product Catalogue

Joint


Spine

Sports Med

INDEX

1. PRODUCT DESCRIPTION	4
2. INDICATION	4
3. REFERENCES AND SIZES	5

1. PRODUCT DESCRIPTION



Cortical Fixation System for tendons and ligaments. The cortical fixation system comprises a continuous loop for graft suspension and a titanium button for the cortical support.

The device is composed of:

- 1 titanium button (ISO 5832-3)
- 1 polyethylene terephthalate loop (PET)
- 2 polyethylene terephthalate pulling and flipping threads (green and white)



Product is delivered in a double peel package and Gamma sterilized.

Extended titanium button available (TA6V ELI). Product is provided in a double peel package and gamma sterilized.

2. INDICATION

The cortical fixation systems are used for fixation of tendons and ligaments during orthopaedic reconstruction procedures such as anterior cruciate ligament (ACL) reconstruction.

3. REFERENCES AND SIZES

REF. NO.	DESCRIPTION	PICTURE
05.05.0080	Continuous Loop 15 mm	
05.05.0081	Continuous Loop 20 mm	
05.05.0082	Continuous Loop 25 mm	
05.05.0083	Continuous Loop 30 mm	
05.05.0084	Continuous Loop 35 mm	
05.05.0085	Continuous Loop 40 mm	
05.05.0086	Continuous Loop Extended Plate	

MectaLoop continuous loop button for ACL reconstruction is a class IIb medical device manufactured by COUSIN BIOTECH s.a.s. The CE conformity evaluation has been carried out by the notified body SGS Belgium NV (CE1639). The quality management system of COUSIN BIOTECH s.a.s. has been certified as meeting the requirements of ISO 13485. Please read instructions for use carefully.

Non contractual pictures and texts. Specifications likely to be modified without notice.

Cousin Biotech S.A.S. capital : 340 656 € - 398 460 261 RCS Lille – N°TVA FR 34 398 460 261



**REDEFINING BETTER
IN ORTHOPAEDICS
AND SPINE SURGERY**

MEDACTA.COM



Distributed by
Medacta International SA
Strada Regina - 6874 Castel San Pietro - Switzerland
Phone +41 91 696 60 60 - Fax +41 91 696 60 66
info@medacta.ch

Find your local dealer at: medacta.com/locations

All trademarks and registered trademarks are the property of their respective owners.
This document is not intended for the US market.
Please verify approval of the devices described in this document with your local Medacta representative.

MectaLoop
Product Catalogue

ref: 99.119SMK.180
rev. 00

Last update: April 2020
CE 1639