

CORTICAL REAMER PIN



DESIGNED TO STREAMLINE ACL RECONSTRUCTION


**CREATE A 4.5 MM TUNNEL
IN A SINGLE STEP**

TECHNICAL DATA

LENGTH: 380 mm
DIAMETER: Ø2.4 mm proximal body &
Ø4.5 mm spade tip

**LASER MARKED SCALE
FOR TUNNEL LENGTH
EVALUATION**

**PROXIMAL EYELET TO USE
THE INSTRUMENT AS A
PASSING WIRE**



ORDERING INFORMATION

Order Number	Description
05.05.10.0150	Ø4.5mm cortical reamer pin

DEDICATED INSTRUMENTS

Order Number	Description	Order Number	Description
05.05.10.0136	Femoral Aimer AM Ø7-8 for cortical reamer pin	05.05.10.0137	Femoral Aimer AM Ø9-10 for cortical reamer pin

RELATED PRODUCTS

Order Number	Description	Order Number	Description
05.05.0091	FairFix Adjustable Button	05.05.0023	MectaFix Continuous Loop 20 mm
05.05.0092	FairFix Adjustable Button Extender	05.05.0024	MectaFix Continuous Loop 25 mm
05.05.0080	MectaLoop Continuous Loop 15 mm	05.05.0025	MectaFix Continuous Loop 30 mm
05.05.0081	MectaLoop Continuous Loop 20 mm	05.05.0026	MectaFix Continuous Loop 35 mm
05.05.0082	MectaLoop Continuous Loop 25 mm	05.05.0027	MectaFix Continuous Loop 40 mm
05.05.0083	MectaLoop Continuous Loop 30 mm	05.05.0028	MectaFix Continuous Loop 45 mm
05.05.0084	MectaLoop Continuous Loop 35 mm	05.05.0029	MectaFix Continuous Loop 50 mm
05.05.0085	MectaLoop Continuous Loop 40 mm	05.05.0030	MectaFix Continuous Loop 55 mm
05.05.0086	Continuous Loop Extended Plate	05.05.0031	MectaFix Continuous Loop 60 mm
05.05.0022	MectaFix Continuous Loop 15 mm		





SPORTSMED.MEDACTA.COM



Distributed by **Medacta International SA** | Strada Regina - 6874 Castel San Pietro - Switzerland | Phone +41 91 696 60 60 | Fax +41 91 696 60 66 | info@medacta.ch

Find your local dealer at: medacta.com/locations. All trademarks and registered trademarks are the property of their respective owners. This document is intended for the US market. Please verify approval of the devices described in this document with your local Medacta representative.

ref: 99.126SMK.1SALESSUS rev. 00 Last update: July 2022