

SnugFit ALL-SUTURE

BRAIDED SOFT ANCHOR



ALL-SUTURE ANCHOR FOR A SECURE AND RELIABLE FIXATION

MULTIPLE ANCHOR CONFIGURATIONS AVAILABLE

**INDICATION: ROTATOR CUFF
& BICEPS TENDON REPAIR**

Double / Triple Loaded USP#2 HS Fiber Suture
Single / Double Loaded USP#2 HS Fiber Suture &
Single Loaded 1.4mm Tape


SMALL ANCHOR DIAMETER

SINGLE, DOUBLE & TRIPLE LOADED

Allowing you to attach more soft tissue

Ø2.8MM DRILL

For minimal bone removal



ORDERING INFORMATION

Order Number	Description	Order Number	Description
05.11.023	SnugFit All-Suture anchor Size 2 - Short - 2x USP#2	05.11.025	SnugFit All-Suture anchor Size 2 - Short - 3x USP#2
05.11.024	SnugFit All-Suture anchor Size 2 - Short - 1x USP#2 / 1x Tape	05.11.026	SnugFit All-Suture anchor Size 2 - Short - 2x USP#2 / 1x Tape

DEDICATED INSTRUMENTS

Order Number	Description	Order Number	Description
05.10.10.0010	Aimer Crown Tip - Short	05.11.10.0063	SnugFit Anchor - Drill Size 2 - Short
05.10.10.0011	Aimer Fishmouth Tip - Short	05.10.10.0013	Obturator Short
05.10.10.0012	Aimer Fork Tip - Short		

RELATED PRODUCTS

Order Number	Description	Order Number	Description
05.15.10.0015	FastShuttle - Suture Passer & Retriever	05.10S.001	Shoulder General Set
05.15.10.0016	FastShuttle - Suture Passer	05.14.10.0001	Shoulder Cannula System Ø8.5x70mm
05.15.10.0017	FastShuttle - Disposable Needle (5ea)	05.14.10.0002	Shoulder Cannula System Ø8.5x90mm
05.15.10.0021	FastShuttle Labral Suture Passer - 45° Left - STERILE (1x)	05.14.10.0035	Shoulder Cannula System Ø8.5x70mm for Ø4mm Switching Stick
05.15.10.0022	FastShuttle Labral Suture Passer - 45° Right - STERILE (1x)	05.14.10.0036	Shoulder Cannula System Ø8.5x90mm for Ø4mm Switching Stick

